Evolve Type-S™ Locating Bar

For use with Evolve

INSTRUCTIONS FOR USE (IFU)

Towards Gantry

Towards Gantry

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Kızığa Doğru

INTENDED USE

The device is indicated to aid in supporting and positioning adult and adolescent patients undergoing radiation therapy of the breast and chest region including electron, photon, and proton treatments. Device is also used to position the patient during image acquisition to support treatment planning including in Computed Tomography (CT) and Magnetic Resonance (MR) imaging systems. The device is not intended for use with patients under 12 years of age.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

**WARNING**

- Do not use if device appears damaged.
- Do not use to secure positioning devices.
- Do not use device with unapproved accessories or devices.
- Do not allow patient to reposition themselves.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Do not use excessive force on device or accessories. Device damage may occur.
- Ensure gantry indicator is positioned towards gantry.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION


- Device is MR safe.

MR Safe

USING EVOLVE™ WITH SLIDE GUIDE 3-PIN LOK-BAR™

1. Locate Lok-Bar™ at appropriate longitudinal position.


**WARNING**


- Ensure Lok-Bar™ is secure.

2. Place Evolve™ onto Lok-Bar™ pins.


NOTE: Ensure device is fully supported.

3. Hook Evolve™ Type-S™ Locating Bar on stabilizer bar.


4. Slide entire system into place.

5. Push T-Pins down to lock Type-S™ Locating Bar in place.

NOTE: Lift T-Pins to unlock.

REPROCESSING**WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, clean surface by removing visual contamination with 70% alcohol, Theracide (PDI Sani-Cloth Plus) wipes, or PDI Super Sani Wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. LOK-BAR AND TYPE-S ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC