ProForm™ HNS Leo Cancer Care Adjustable Backrest

For use with Leo Cancer Care (MARIE™, RUBY™, AIDA™)

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

- The device is used to assist with immobilization, positioning and reposition patients undergoing radiation therapy.
- The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Do not apply a load exceeding 363 lbs (165 kg) to the entire system including accessories.
- Do not apply a load exceeding 51 lbs (23.3 kg) to backrest.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at www.CQmedical.com for WET values.
- Do not allow patient to reposition themselves.
- Do not reposition device with patient on it.
- Do not use excessive force on device. Injury or device damage may occur.


NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

INSTALLING BACKREST**⚠ WARNING**

- Always verify device characteristics prior to treating patients.


1. With both hands, grab each side of backrest (1) and stand in front of seat pan (7).
2. Orient backrest so Leo Cancer Care logo faces the patient side.
3. Line mounting pin (3) up with backrest interface (4) so mounting pin will insert into mounting holes (5).
4. Lower backrest until mounting pin is firmly seated in mounting hole.

- NOTE:** Ensure backrest is in docked position by pushing downward equally on top of backrest.
5. When fully seated, the control system allows continued function of the system. If it does not, refer back to step 4.

**ADJUSTING EXTENSION**

1. To adjust extension: Depress handles (6) on both sides to slide extension (2) up and down (Figure 1).
2. Release handles when desired height is reached.
3. Apply slight pressure up and down to ensure handle engagement is locked in place.
4. Once height adjustment is complete, lock handles by pulling up (away from the backrest) to prevent movement (Figure 2).


5. Position patient as dictated by setup protocol.

REMOVING BACKREST

1. Stand in front of seat pan.
2. With both hands, grab each side of backrest (1) and lift upward until backrest is removed from mounting holes (2).

NOTE: A slight resistance will be felt as the ball catch releases from mounting hole.

3. Set aside backrest.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with 70% alcohol or Theracide (PDI Sani-Cloth Plus) wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. PROFORM IS A TRADEMARK OF MEDTEC LLC.
ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**