ONEClamp

For use with Body Pro-Lok™ Respiratory Belt and ZiFix Traverse™ Abdominal/Thoracic Motion Control System

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
Kizağa Doğru

INTENDED USE

- The device is indicated to assist in the proper positioning of patients for radiation therapy including electron, photon and proton treatments, and simulation including CT and MR image acquisition.
- The device is intended to apply abdominal compression for managing movement and tumor motion during respiration while maintaining maximum comfort to the patient.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Inspect device Velcro for damage, wear or contamination that may cause belt to loosen over time.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Ensure patient is set up in correct position for each treatment before beginning treatment.
- Ensure patient remains still for duration of setup and treatment.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between devices and all linac or imaging equipment before use.
- Remove sharp objects from patient and treatment area to avoid puncturing.
- Pressure indication is for reference only.
- Clinical personnel must monitor patient constantly when using respiration suppression accessories to ensure patient can breathe adequately.
- Ensure device is not too tight for patient comfort.
- Standing or sitting on ZiFix Traverse™ Abdominal/Thoracic Motion Control System may damage clamp or cause injury.
- Do not carry device using ZiFix Traverse™ Abdominal/Thoracic Motion Control System.
- Ensure clamps are properly latched.
- Use only CQ Medical supplied hand pump.

- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - The ZiFix Traverse™ Abdominal/Thoracic Motion Control System will attenuate a radiotherapy beam. Actual attenuation based on setup should be verified with your particular equipment. Attenuation and increased skin dose should be taken into account during planning and treatment.

MRI SAFETY INFORMATION

Body Pro-Lok™ Respiratory Belt and ONEClamp

- Non-clinical testing has demonstrated the Respiratory Belt is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.5°C after 15 minutes of continuous scanning.
- Accessory Hand Pump with Pressure Gauge is MR Conditional. Refer to manufacturer's instructions for conditions of safe use.


ZiFix Traverse™ Abdominal/Thoracic Motion Control System

- Non-clinical testing has demonstrated the Manometer Pump is MR Conditional. This device may be used in an MR system under the following conditions:
 - Static magnetic field of 3T or less.
 - Disconnect the Manometer Pump from the device before scanning and remove the Manometer Pump from the patient table.
 - Prior to scanning, move the manometer pump into gradient fields no greater than 200 mT (2000G), i.e. weaker gradient sections, of the magnetic field of the MRI scanner.
 - If the magnetic field strength lines in the room configuration are not known, disconnect and remove device from the MR room prior to scanning. Consult the user manual for the MR system for information regarding magnetic field strength and spatial gradients for the system.
- Non-clinical testing has demonstrated the ZiFix Traverse Abdominal/Thoracic Motion Control System is MR Conditional. Once the conditions for the Manometer Pump are met, the ZiFix Traverse Abdominal/Thoracic Motion Control System may be used in an MR system meeting the following conditions:
 - Static magnetic field of 3T or less.


ADJUSTING ONECLAMP

1. Insert a thin, sturdy object such as a paperclip, into pinhole and raise or lower adjustable width clips as needed. Set to Configuration 1 for 4 inch belts and Configuration 2 for 5 inch belts.


NOTE: In the event ONEClamp stops responding to knob, hold release button while turning knob.

USING BODY PRO-LOK™ RESPIRATORY BELT

NOTE: Indicator bar can be placed in superior or inferior position. Ensure indicator bar orientation is consistent between use of all belt systems.

1. Attach clamps to side rails. Secure by tightening knob until it clicks.
2. Place Velcro straps loosely through each clamp.

NOTE: Ensure adjustable width clips of clamp are set for correct belt width.


⚠ WARNING

- Ensure adjustable width clips are engaged.

3. Hold respiratory belt firmly in place over desired area.
4. Place Velcro strap over respiratory belt.
5. Place second Velcro strap over respiratory belt, covering first strap if necessary.
6. Attach respiratory gauge and pump to respiratory belt.
7. On gauge, use range A through H for respiratory belt inflation.
8. Inflate respiratory belt until desired respiratory restriction is achieved.
9. Prior to imaging or treatment, remove gauge and pump from respiratory belt.

NOTE: Note all indexed positions, including orientation of indicator bar, and acceptable pressure range on setup sheet.

USING ZIFIX TRAVERSE™ ABDOMINAL/THORACIC MOTION CONTROL SYSTEM


INITIAL ASSEMBLY AND POSITIONING ON PATIENT

1. Ensure paddle is appropriately aligned for attachment to the bladder. Align hook-and-loop fastener.


NOTE: • If using white alignment lines on paddle, complete alignment prior to attaching bladder. Alignment lines are for reference use only.
 • The virtual indexing marks may be used to verify location in imaging. The virtual indexing marks are for reference use only.

2. Position belt with numbers facing down. Feed the long end of belt through clamp. Return belt on itself ensuring full length of hook is engaged on loop.

NOTE: Ensure adjustable width clips of clamp are set for correct belt width.

- Position paddle and bladder on patient with tip of paddle aligned with patient's sternum and below xiphoid process.

NOTE: The five (5) virtual indexing marks located around paddle perimeter may be used to aid in patient setup relative to patient anatomy under CT or X-ray imaging.


- Attach clamp. Secure by tightening knob until it clicks.

NOTE: Ensure there is no material between the clamp and device when installing.

- Lay belt over patient, ensuring number side is facing down and gantry indicator points towards gantry, attaching it to paddle and bladder. Ensure loop on belt aligns and engages with hook on bladder.
- Install second clamp on opposite device edge. Secure by tightening knob until it clicks.
- Tighten belt as patient tolerance will allow.
- Record clamp and belt scale positions.

INFLATING BLADDER

- Close valve by rotating red knob towards you. Ensure shutoff valve is open by checking valve handle is in-line with tubing.


- Inflate bladder to desired pressure by squeezing bulb repeatedly.

⚠ WARNING

- Do not inflate bladder pressure beyond 125 mmHg.

- Close shutoff valve by turning valve handle perpendicular to tubing. Wait 2 minutes.

NOTE:

- It is recommended to reinflate bladder to initial pressure after 2 minutes of initial inflation by opening the shutoff valve, squeezing the manometer bulb, and shutting the shutoff valve.
- During normal operation, fluctuation of 10 mmHg, can be expected.

- After 2 minutes, open shutoff valve by turning valve handle in-line with tubing.
- Inflate bladder to desired pressure by squeezing bulb repeatedly. Close shutoff valve.

NOTE:


- Record manometer readout.
- This device is to be used as reference only and should not be used in any other application.

⚠ WARNING

- Max bladder pressure is 125 mmHg.
- Pump must be disengaged via "Quick Disconnect Fitting" and removed from patient table before scanning or treatment.

REMOVING PUMP FROM BLADDER

- Verify shutoff valve is perpendicular to tubing.
- Disconnect pump by pressing the grey button on Quick Disconnect Fitting.

NOTE: Expect a small decrease in pressure if pump is reconnected and valve is reopened.

REMOVING BELT FROM TABLE

- Disconnect Quick Disconnect Fitting. Open shutoff valve and release the air from bladder. OR rotate red knob away from you.
- Loosen belt as needed. Unlock clamp and firmly slide off edge.

REPROCESSING **WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- User must clean device between uses according to instruction.

Body Pro-Lok™ Respiratory Belt

1. If necessary, clean surface by removing visual contamination with 70% alcohol or Theracide (PDI Sani-Cloth Plus) wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

NOTE: Color may transfer to wipe during normal cleaning.

ZiFix Traverse™ Abdominal/Thoracic Motion Control System

The ZiFix Traverse™ Abdominal/Thoracic Motion Control System can be cleaned using standard cleaning practices and common WATER-BASED hospital cleaning agents. The following cleaning materials have been tested and found to be appropriate:

- Cidex® 2.4% Activated Dialdehyde Disinfecting Solution
- Soap and Water
- Isopropyl Alcohol

NOTE: • DO NOT use other solvents to clean the device.
• The gauge should be cleaned with a soft, dry cloth.

MAINTENANCE **WARNING**

- Inspect device prior to use for signs of damage and general wear.

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CEpartner4U
Esdoornlaan 13
3951 DB Maarn, Netherlands
Ph: +31 (0) 6-516.536.26

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

 MEDTEC LLC
1401 8th Street SE
Orange City, IA 51041
United States
info@cqmedical.com

COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO LOK AND ONECLAMP ARE TRADEMARKS OF MEDTEC LLC. INFINITYEDGE, KVUE, LOW-E, AND ZIFIX TRAVERSE ARE TRADEMARKS OF QFIX. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.

www.CQmedical.com