Posi Indexing Bar

For use with MTM3303, MTM3304, MTM3305 and MTM3004

INSTRUCTIONS FOR USE (IFU)

Towards Gantry

Towards Gantry

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Kızığa Doğru

INTENDED USE

The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

 WARNING

- Do not use if device appears damaged.
- Use MR compatible accessories only.
- Only use CQ Medical compatible accessories.
- Ensure all accessories used during simulation are used during treatment.
- Ensure accessories are securely attached in correct position prior to imaging.
- Ensure device is secure prior to use.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Indexing Bar is a positioning device only, do not use to secure patient.
- Do not allow patient to reposition themselves.
- Verify patient position with completed setup sheet prior to treatment.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

MR Safe

- Device is MR safe.
- Image artifacts may be present within the imaging region; evaluate and adjust as necessary.


USING INDEXING BAR


 WARNING

- Ensure gantry indicator is positioned towards gantry.

NOTE: Device is compatible with CQ Medical - MTM3303, MTM3304, MTM3305, MTM3004.

STEP 1

STEP 2**REPROCESSING****WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.