
Encompass™ Shoulder Thermoplastic

INSTRUCTIONS FOR USE (IFU)

CE

PARTS LIST	
RT-B889KYSH	Encompass™ Fibreplast® Shoulder Thermoplastic
RT-B889KYSH-1	Encompass™ Shoulder Thermoplastic Adapter

INTENDED USE

The Encompass™ SRS Immobilization System provides noninvasive stereotactic head and neck immobilization by using a patient specific thermoplastic mask that conforms to the patient's features to provide accurate, reproducible positioning, repositioning, and immobilization. The Encompass™ SRS Immobilization System allows the patient to undergo diagnostic imaging in the same position as that of treatment, enabling radiation therapy.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- No modification of this equipment is allowed. If any part of this device experiences a catastrophic load, appears damaged or functions improperly, discontinue use immediately and contact CQ Medical at info@CQmedical.com.
- To ensure safe use of device, it is recommended that the user is educated and trained on safe operation of system prior to use.
- Only use CQ Medical compatible accessories.
- Thermoplastic material may cause skin irritation in individuals with sensitive skin.
- Cutting thermoplastic masks may reduce the rigidity of the device. If cutting is necessary, ensure rigidity of mask is maintained by adding thermoplastic strips. Thermoplastic strips are not included with masks but can be purchased by contacting your CQ Medical representative.
- Do not sterilize thermoplastic.
- Thermoplastic masks will attenuate a radiotherapy beam and increase skin dose. Attenuation and increased skin dose should be taken into account during planning and treatment.
- Dose depth, deposition, and transition area effects must be evaluated during planning and treatment within a proton therapy environment.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at www.CQmedical.com for WET values.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.


- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Refer to www.CQmedical.com for a listing of symbols and their definitions.

MRI SAFETY INFORMATION

	Non-clinical testing has demonstrated the following components are MR Safe. • Thermoplastic Mask (RT-B889KYSH) • MR image quality may be compromised if area of interest is relatively close to edge of mask. Optimize MR imaging parameters prior to treating patient.
	Non-clinical testing has demonstrated the following components are MR Conditional. These devices may be used with an MR system meeting the following conditions: • Shoulder Adapters (RT-B889KYSH-1) • Static magnetic field of 1.5 T and 3.0 T • Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)


USING SHOULDER ADAPTERS

1. Install Shoulder Adapter as desired on top surface of Encompass™. Shoulder Adapters can be installed in either a wide or narrow orientation.

2. Install both barrel nuts with included brass screw and secure Shoulder Adapter to Encompass™.

NOTE: Ensure device is securely attached.

**USING THERMOPLASTICS****⚠ WARNING**

- Ensure thermoplastic is properly aligned over patient before forming.
- Ensure patient is comfortable during molding process to ensure fit is comfortable during future use.
- Form thermoplastic to avoid contact with tracheotomy tube.
- When securing mask, take care to not pinch patient or operator between mask and device.
- Ensure thermoplastic is properly attached using pins.

Thermoplastic Type	Waterbath Temperature	Waterbath Heating Time	RapidHeat™ Oven Temperature	RapidHeat™ Oven Heating Time
Encompass™ Fibreplast® Thermoplastic	74°C (170°F)	2 minutes	74°C (165°F)	7-8 minutes

1. Setup waterbath with 2–3" of water and heat or turn on dry heat oven (refer to table for appropriate temperature).
2. Place thermoplastic into waterbath/dry heat oven for directed time.

NOTE: Improper heating time may cause difficult application.

3. Remove thermoplastic from waterbath/dry heat oven.
4. Quickly remove excess water from thermoplastic (*less than five seconds*).


WARNING

- Water is hot, do not reach into water to remove device.
- Feel thermoplastic mask to make sure it is not too hot to come in contact with patient's skin.


5. Position thermoplastic over treatment area.
6. Gently pre-stretch thermoplastic.
7. Gently pull sides of frame outward and down to attachment base.

NOTE:

- Some thermoplastic masks may require assistance in pulling of mask.
- Overstretching of thermoplastic material may affect rigidity of mask. After pulling thermoplastic mask, ensure rigidity of mask is maintained.
- While pulling mask ensure thermoplastic is not folded incorrectly under itself.

8. Secure thermoplastic by pressing down on all available Shoulder Frame pins. Ensure all pins are engaged in adapter.

9. Gently mold material to contours of treatment area.
 10. Allow thermoplastic to cool completely before releasing from attachment base.
- NOTE:
- Removing thermoplastic before completely cool may cause thermoplastic deformations, such as shrinkage or collapse.
 - Use of an immobilization device which aids shoulder depression, such as patient handles, is recommended when using Shoulder Thermoplastic.

REPROCESSING


WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface of thermoplastic and/or frames and pins by removing gross visual contaminants with alcohol or warm soapy water. If gross visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2026 ALL RIGHTS RESERVED. CQ MEDICAL IS A REGISTERED TRADEMARK OF MEDTEC LLC. ENCOMPASS IS A TRADEMARK OF QFIX. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC