Upright Arm Positioner

For use with Leo Cancer Care Systems and Mevion S250-FIT Proton Therapy System™

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is designed as a support and fixation device for patients undergoing radiotherapy treatment(s).

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING


- No modification of this equipment is allowed. If any part of this device experiences a catastrophic load, appears damaged or functions improperly, discontinue use immediately and contact CQ Medical at support@CQmedical.com.
- Ensure device is securely attached.
- Ensure accessories are securely attached in correct position prior to imaging.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Do not apply a load exceeding 66 lbs (30 kg).
- Do not allow patient to reposition themselves.
- Do not reposition device with patient on it.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


INSTALLING UPRIGHT ARM POSITIONER

NOTE: Backrest (consult appropriate Instructions for Use) must be installed before installing upright arm positioner.


1. Position patient with Respiratory Belt and other desired positioning accessories (consult appropriate Instructions for Use) per patient needs.
2. Determine proper orientation of right hand and left hand clamps as shown.


3. Install clamp by sliding onto backrest from side. If necessary, loosen clamp knob by rotating counterclockwise to allow clamp to fully seat on backrest rail and indexing pin.


4. With one hand, hold clamp at desired location against backrest rail. With other hand, tighten clamp knob by rotating clockwise until clamp audibly clicks.


5. Repeat Step 4 for remaining side. Ensure clamps are aligned at matching index positions. Record index position of clamps.

6. Locate clamp slots on both clamps. Install upright arm positioner by inserting arm positioner rails at a downward angle through clamp slots. Drop down to set position. Record index position.


7. Place posts in position on upright arm positioner per patient needs. Record index position of both posts.

REMOVING UPRIGHT ARM POSITIONER

1. Remove posts from upright arm positioner.
2. Lift upright arm positioner at an angle and slide rails from clamps.
3. Loosen clamp knobs until clamps can be removed from either side. In the event the clamp stops responding to knob, hold release button while turning knob.

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with 70% alcohol or Theracide (PDI Sani-Cloth Plus) wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2026 ALL RIGHTS RESERVED. CQ MEDICAL IS A REGISTERED TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**