AirShuttle™ Transfer Handles

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

The Symphony® AirShuttle™ Transfer Handles are indicated to aid in the support, positioning, and transfer of a patient for procedures involving imaging, including MRI; and external beam radiation therapy treatment with electrons, photons or protons; and other procedures requiring transfer of a patient. The Symphony® AirShuttle™ Transfer Handles are designed to interface with Symphony® AirShuttles.

PATIENT TARGET GROUPS

Patients undergoing radiation therapy, diagnostic imaging procedures, or other procedures involving transfer of a patient.

INTENDED USERS

Person qualified in accordance with requirements of the regulatory region.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- No modification of this equipment is allowed. If any part of this device experiences a catastrophic load, appears damaged or functions improperly, discontinue use immediately and contact CQ Medical at support@CQmedical.com.
- Patient should be examined for any pre-existing conditions that might prohibit safe use of device.
- To prevent patient injury, ensure device is positioned for proper patient positioning in accordance with facility, procedure, and provider requirements.
- Ensure all accessories including the transfer handles are secured to the AirShuttle™ prior to use.
- The AirShuttle™ Transfer Handles must be removed from the AirShuttle™ prior to placing the AirShuttle™ into the bore of any imaging system.

NOTE:

- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
- Refer to www.CQmedical.com for a listing of symbols and their definitions.

MRI SAFETY INFORMATION

Use only proven MR Safe or MR Conditional accessories, tested, and approved for your MR system. Consider MR compatibility of accessories prior to use on MR system.

⚠ WARNING

- Use of unapproved MR accessories may result in:
 - Injury to Patient
 - Patient burns
 - Damage to equipment
- Recommended maintenance and service, as well as use of only CQ Medical-provided accessories, components, and replacement parts, are necessary to ensure safety, performance, and MRI compatibility of product(s), as well as to maintain applicable warranties.


MR Conditional

Non-clinical testing has demonstrated this device is MR Conditional and may be used with an MR system meeting the following conditions:


- Static magnetic field of 3T or less.

INSTALLATION**ATTACHING THE TRANSFER HANDLES**

1. Locate Transfer Handle mounting points on the desired AirShuttle™.
2. Slide locking component of the Transfer Handles into mounting points on the AirShuttle™.
3. Once installed, pull on Transfer Handles to ensure they are properly locked.

**REMOVING THE TRANSFER HANDLES**


1. Squeeze release trigger on either side of Transfer Handles, then pull handles from AirShuttle™ mounting points.

MAINTENANCE

NOTE: Do not submerge device in liquid.

- 1. Device can be cleaned with a mild, non-abrasive cleaning or disinfecting solution. To clean, apply solution to clean cloth and wipe surface. Visually inspect device, if it is not clean, repeat previous cleaning steps until visually clean. Use a clean cloth moistened with water to wipe device to remove any cleaning agent residue. To dry, wipe device with a clean, dry cloth. The following cleaning material has been tested and found to be appropriate for cleaning device.
 - Water
 - Soap and Water
- 2. The following products are recommended. Refer to specific instructions from chemical agent manufacturer.
 - 10% Clorox® Bleach Solution
 - Isopropyl Alcohol

SPECIFICATIONS	
MECHANICAL SPECIFICATIONS	
Overall Dimensions	RT-5100-10-B: 644mm x 200mm x 152mm (25.4in x 7.9in x 6in) RT-5100-10: 441mm x 200mm x 152mm (17.4in x 7.9in x 6in)
Overall Weight	RT-5100-10-B: 1.43kg (3.15lbs) RT-5100-10: 1.27kg (2.8lbs)

NOTE: Lithotomy version of Transfer Handles shown.

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
ÖZELLİKLER	
MEKANİK ÖZELLİKLER	
Genel Boyutlar	RT-5100-10-B: 644mm x 200mm x 152mm (25,4inç x 7,9inç x 6inç) RT-5100-10: 441mm x 200mm x 152mm (17,4inç x 7,9inç x 6inç)
Genel Ağırlık	RT-5100-10-B: 1,43kg (3,15lbs) RT-5100-10: 1,27kg (2,8lbs)

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 QFix
440 Church Rd
Avondale, Pennsylvania
United States
info@cqmedical.com

COPYRIGHT © 2025 ALL RIGHTS RESERVED. CQ MEDICAL IS A REGISTERED TRADEMARK OF MEDTEC LLC. AIRSHUTTLE IS A TRADEMARK OF QFIX. SYMPHONY IS A REGISTERED TRADEMARK OF QFIX. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.

www.CQmedical.com