PRODUCT GUIDE & USER MANUAL

RT-4550KVGE, RT-4550KVGEL, RT-4550KVGER,
RT-4550KVPBR, RT-4550KVSS1, RT-4550KVSS2,
RT-4550KVSS3, RT-4550KVSS4, RT-4550KVSS5,
RT-4550KVTAQ, RT-4550KVTAQ2, & RT-4450KVPH01

kVue™ CT Overlay

 EC REP Advena Limited Tower Business Centre, 2nd Flr Tower Street, Swatar, BKR 4013 Malta	UK Responsible Person QServe Group UK, Ltd. 49 Greek Street Soho, London W1D 4EG United Kingdom
 Made in the USA by Qfix 440 Church Rd, Avondale, Pennsylvania, USA +1 484-720-6053 www.Qfix.com	CH REP Raditec Medical AG Schlossberg 5a 5454 Bellikon Switzerland

2002600_ML_1 / 2024-08

EN	English	KK	Kazakh (Қазақша)
DA	Danish (Dansk)	NL	Dutch (Nederlands)
DE	German (Deutsch)	RO	Romanian (Română)
ES	Spanish (Español)	RU	Russian (Русский)
FR	French (Français)	SV	Swedish (Svenska)
IT	Italian (Italiano)	ZH	Chinese (简体中文)
JA	Japanese (日本語版)		

EN: English

GENERAL PRECAUTIONS

WARNING STATEMENTS

! WARNING ! NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED. IF ANY PART OF THIS DEVICE EXPERIENCES A CATASTROPHIC LOAD, APPEARS DAMAGED OR FUNCTIONS IMPROPERLY, DISCONTINUE USE IMMEDIATELY AND CONTACT QFIX AT +1 484-720-6053 OR TECHSUPPORT@QFIX.COM.

! WARNING ! THE ONETOUCH™ IS A PRECISION MECHANISM AND IS ONLY FOR USE WITH APPROVED QFIX KVue DEVICES. UNAPPROVED NON-QFIX DEVICES THAT HAVE NOT BEEN VALIDATED BY QFIX MAY NOT BE SAFE FOR USE AND WILL VOID WARRANTY IF USED.

SERIOUS INCIDENTS

Please report any serious incidents (e.g. incidents which result in or have the potential to result in death or serious injury) to both Qfix and your country's Competent Authority.

ADDITIONAL WARNINGS

- The kVue CT will not work with certain panels with all inserts, such as those with a beam limiter.
- Pinch points may exist and are properly marked on the kVue CT Overlay.
- Attention is necessary when moving the Couch and CT Insert into the Gantry Bore to ensure that pinch point risk to hands and fingers are avoided. If cradle is retracted from the bore of the machine with the cradle in its lowest position, it may be possible to damage the front end eccentric cams.
- The unit must be inspected periodically for damage and atypical wear.
- Use of this product in a manner not specified in the User Manual may be unsafe and is not recommended.
- Use of unapproved devices or attachments may be unsafe and is not recommended.

! CAUTION ! LOCATING BARS USED IN CT SCANS MAY CAUSE ARTIFACTING. MOVE LOCATING BAR OUTSIDE REGION OF INTEREST IF NECESSARY.

LOAD RATING

DO NOT exceed 249 kg (550 lbs) or the lowest of the kVue CT Overlay, kVue, and CT weight limit.

! NOTE ! When combined with the OEM couch base, the load rating is the lower of the two safe working loads. The kVue CT Overlay load should NOT exceed the original couch base manufacturer's specifications. Refer to the product literature provided by the original manufacturer.

WARNING LABELS & DESCRIPTIONS

Refer to Qfix.com for a listing of symbols and their definitions.

	LOAD RATING DO NOT exceed 249 kg (550 lb) uniformly distributed load.
	PINCH POINTS Pinch points may exist and are indicated by one of these symbols.

INTENDED USE

This device is intended to immobilize, position and reposition patients undergoing radiation therapy.

! NOTE ! United States Federal law restricts this device to sale by or on the order of a physician.

PATIENT TARGET GROUPS

Patients undergoing radiation therapy or diagnostic imaging procedures.

INTENDED USERS

The intended user for the products is a person qualified in accordance with the requirements of the regulatory region.

FEATURES

DESCRIPTION

The kVue CT is a CT simulator that allows kVue Inserts to be placed directly into the CT. The dosimetric properties of your couch top can now be taken into account in the Treatment Planning Software (TPS). The kVue CT is available in many CT makes and models.

RT GATING CAMERA

Location for mounting RT Gating Camera (Fig. 1).


Fig. 1

TREATMENT INFORMATION

KVUE INSERT

kVue Inserts are constructed from low-density composite materials and will attenuate the photon beam by approximately 1% for AP/PA treatments at 6 MV. Actual attenuation based on setup should be verified with your particular equipment. Treatment through any device, even one constructed on composite materials, will lead to an increase in skin dose.

OPERATING INSTRUCTIONS

INSTALLATION

KVUE COUCH TOP INSTALLATION

The kVue Couch Top is designed to fit all major radiotherapy treatment systems. Quick installation minimizes linear accelerator downtime during installation. Refer to the kVue Couch Top Installation Manual for specific installation information.


! NOTE ! As detailed in the kVue Couch Top Installation Manual, the final height of the kVue surface may be different from the original system. Refer to the original system instructions to reset the system zero (origin) prior to initial use.

KVUE INSERT INSTALLATION

1. Place the kVue Insert on top of the support beams and align the two mating pins with the receiving holes in the OneTouch Latch.
2. From the head end of the kVue, slide the kVue Insert directly into the receiving holes.

You will hear a click when the kVue Insert is locked into place. **ONLY** the green button surface will be visible when properly installed.

3. Verify insert is fully engaged by gently pulling it backwards, away from OneTouch Latch, to ensure it does not come out of receiving holes. Repeat this step during patient setup, even if insert has previously been installed.

KVUE CT INSTALLATION

! NOTE ! The following photos show installation of **QUANTUM™ CT**.

1. Place the kVue CT Overlay on the top of the CT Couch centered left to right and overhanging the CT Head by 5 to 10 cm (2 to 4 inches) (Fig. 2).


Fig. 2

2. Slide the kVue CT Overlay towards the CT Foot End keeping the head and foot ends near to the centerline of the CT couch (Fig. 3).


Fig. 3

3. The handle will move to the latched position automatically when the insert has been engaged properly. When only the green shows and no red, the kVue CT Overlay is locked on (Fig. 4).


Fig. 4

SET-UP

LOCATING BAR INSTALLATION

The Locating Bar has two locating pins that match most standard positioning accessories. To attach, place either end of the Locating Bar in the appropriate Varian Exact® compatible index notches and snap in place.

REMOVAL

LOCATING BAR REMOVAL


Pull up on either end of the Locating Bar.

KVUE INSERT REMOVAL

! WARNING ! NEVER LIFT THE END OF THE KVUE INSERT ! ALTHOUGH THE MECHANISM IS VERY ROBUST, THE MATING PINS MAY BEND OR BREAK, RENDERING THE KVUE UNUSABLE.

1. Rotate Lever counter-clockwise.
2. Pull the rotated Lever from one side of the kVue until the kVue Insert is released.
3. Slide the kVue Insert away from the kVue.

! NOTE ! The images shown below are of a standard kVue Couch Top. kVue Insert removal is the same for all versions of the kVue.

REMOVAL

KVUE CT REMOVAL

1. Pull/push the release handle on the KVue CT Overlay, in the direction indicated by the unlock symbols, to disengage the Overlay from the CT couch (Fig. 5).

Fig. 5

2. While activating the handle, slide the KVue CT Overlay towards the head end 5 to 10 cm (2 to 4 inches) to clear the head end attachments (Fig. 6).

Fig. 6

3. Lift off the KVue CT Overlay.

MAINTENANCE

CLEANING THE SYSTEM

The device can be cleaned with a mild, non-abrasive cleaning or disinfecting solution. To clean, apply solution to clean cloth and wipe the surface. Visually inspect the device, if it is not clean, repeat the previous cleaning steps until visually clean. Use a clean cloth moistened with water to wipe the device to remove any cleaning agent residue. To dry, wipe the device with a clean, dry cloth. The following cleaning material has been tested and found to be appropriate for cleaning the device.

- Water
- Soap and Water

DISINFECTING THE SYSTEM

To disinfect the device surface, refer to specific instructions from the cleaning agent manufacturer.

- 10% Clorox® Bleach Solution
- Isopropyl Alcohol
- Cidex® 2.4% Activated Dialdehyde Solution

DO NOT place sharp objects on the kVue CT.

Periodically check all fasteners for tightness.

CLEANING AND LUBRICATING PROCEDURE FOR ONETOUCH LATCH

We recommend performing this procedure monthly.

Materials List: fluoropolymer (PTFE) dry film lubricant, cotton swabs, and dry wipes or paper towels

! NOTE ! The images shown below are of a standard kVue Couch Top. Maintenance is the same for every kVue Couch Top.

1. Remove kVue Insert (Fig. 7).
2. Wipe any lint or debris from the outside of latch (Fig. 8).
3. Use dry cotton swab to clean out the inside of the receiving holes (Fig. 9).
4. Spray lubricant into the receiving holes (Fig. 10).
5. Use dry cotton swab to wipe out loose lint and debris (Fig. 9).
6. Spray lubricant into the receiving holes again (Fig. 10).
7. Spray lubricant into the gap between the top of the button and the latch housing to lubricate interior components (Fig. 11).
8. Reinsert the kVue Insert and eject by pulling latch from opposite side. Repeat lubrication step 7.
9. Cycle kVue Insert installation and removal 3 to 4 times and repeat lubrication process if Insert does not eject freely.
10. Wipe off any excess lubricant from the exterior of the OneTouch Latch.


Fig. 7


Fig. 8


Fig. 9


Fig. 10


Fig. 11

CLEANING AND LUBRICATING PROCEDURE FOR CT LATCH

We recommend performing this procedure every 3 months.

Materials List: fluoropolymer (PTFE) dry film lubricant, cotton swabs

1. Use dry cotton swab to clean out the inside of the receiving hole (Fig. 12). The hole on the KVue CT is located on the underside of the CT Latch (Fig. 13).
2. Spray lubricant into the receiving holes. (Fig. 14). The hole on the KVue CT is located on the underside of the CT latch (Fig. 15).
3. Use dry cotton swab to wipe out loose lint and debris (Fig. 12 & 13).
4. Spray lubricant into the receiving holes again (Fig. 14 & 15).

Fig. 12


Fig. 13

Fig. 14


Fig. 15

PARTS LIST

- 1. 2000382 – Tall Eccentric
- 2. 2000707 – Short Eccentric (RT-4550SS2 only)
- 3. 2000384 – Guide Block
- 4. 3000118 – Black FHS 3/8-16 x 1"

DEVICE LIST

Part Number	Device Description
RT-4550KVGECE	kVue CT Base for GE Lightspeed VCT Scanner
RT-4550KVGEL	kVue CT Base for GE Lightspeed and HiSpeed CTi CT
RT-4550KVGER	kVue CT Base for GE Lightspeed RT CT
RT-4550KVPBR	KVue CT Base for Philips Brilliance CT
RT-4550KVSS1	kVue CT for Siemens Somatom Plus 4, Volume
RT-4550KVSS2	kVue CT for Siemens Somatom Balance, Emotion
RT-4550KVSS3	kVue CT Base for Siemens Somatom Definition
RT-4550KVSS4	kVue CT Overlay, Siemens
RT-4550KVSS5	kVue CT Overlay, Siemens Biograph PET/CT
RT-4550KVTAQ	KVue CT Base for Cannon Aquilion Scanner for 215 cm cradle
RT-4550KVTAQ2	KVue CT Base for Cannon Aquilion Scanner for 244 cm cradle
RT-4550KVPH01	KVue CT Base for Philips Spectral CT 7500 Scanner

kVue, OneTouch, and QUANTUM are trademarks of Qfix.
Aquilion is a registered trademark of Canon Medical Systems Corporation.
Exact is a registered trademark of Varian Medical Systems.
Cidex is a registered trademark of Johnson & Johnson.
Clorox is a registered trademark of The Clorox Company.

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RT-4550KVGECE, RT-4550KVGEL, RT-4550KVGER,
RT-4550KVPBR, RT-4550KVSS1, RT-4550KVSS2,
RT-4550KVSS3, RT-4550KVSS4, RT-4550KVSS5,
RT-4550KVTAQ, RT-4550KVTAQ2 og RT-4450KVPH01

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