PRODUCT GUIDE & USER MANUAL

RT-4525TILT Tilt-Pro™ Base
RT-4525TILTMRI Tilt-Pro™ MRI Base

 EC REP Advana Limited Tower Business Centre, 2nd Flr Tower Street, Swatar, BKR 4013 Malta	UK Responsible Person QServe Group UK, Ltd. 49 Greek Street Soho, London W1D 4EG United Kingdom
 <i>Made in the USA by</i> Qfix 440 Church Rd, Avondale, Pennsylvania, USA +1 484-720-6053 www.Qfix.com	CH REP Raditec Medical AG Schlossberg 5a 5454 Bellikon Switzerland

2003250_ML_N / 2024-01

EN	English	NL	Dutch (Nederlands)
DE	German (Deutsch)	RO	Romanian (Română)
ES	Spanish (Español)	RU	Russian (Русский)
FR	French (Français)	SK	Slovak (Slovenský)
HU	Hungarian (Magyar)	SV	Swedish (Svenska)
IT	Italian (Italiano)	TR	Turkish (Türkçe)
JA	Japanese (日本語)	ZH	Chinese (简体中文)
KK	Kazakh (Қазақша)		

EN: English

GENERAL PRECAUTIONS

WARNING STATEMENTS

! WARNING ! NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED. IF ANY PART OF THIS DEVICE EXPERIENCES A CATASTROPHIC LOAD, APPEARS DAMAGED OR FUNCTIONS IMPROPERLY, DISCONTINUE USE IMMEDIATELY AND CONTACT QFIX AT +1 484-720-6053 OR TECHSUPPORT@QFIX.COM.

ADDITIONAL WARNINGS

- Always secure the Tilt-Pro to the treatment table prior to use.
- Do not modify the Tilt-Pro in any way.
- Do not use the Tilt-Pro if it appears to have been modified or damaged in any way.
- Never put a patient's full weight on the device. The Tilt-Pro is designed to support only the weight of a patient's head.
- Do not allow the patient to lean on this device.

SERIOUS INCIDENTS

Please report any serious incidents (e.g. incidents which result in or have the potential to result in death or serious injury) to both Qfix and your country's Competent Authority.

MRI SAFETY INFORMATION

 Non-clinical testing has demonstrated the **Tilt-Pro MRI Base, RT-4525TILTMRI**, is MR Safe. This device may be used in an MR system.

TREATMENT BEAM ATTENUATION

The device will attenuate the photon beam. Treatment through the couch top will increase the overall attenuation. Treatment through any device, even one constructed of carbon fiber, will lead to an increased skin dose. Actual attenuation based on setup should be verified with your particular equipment.

WARNING LABELS & DESCRIPTIONS

Refer to Qfix.com for a listing of symbols and their definitions.

MR SAFE

When present, item poses no known hazard in all MR environments. DO NOT use with any unapproved components or accessories. DO NOT use in conjunction with any ferromagnetic or conductive items.

! NOTE ! RT-4525TILTMRI only

INTENDED USE

This device is intended to immobilize, position and reposition patients undergoing radiation therapy.

! NOTE ! United States Federal law restricts this device to sale by or on the order of a physician.

PATIENT TARGET GROUPS

Patients undergoing radiation therapy or diagnostic imaging procedures.

INTENDED USERS

The intended user for the products is a person qualified in accordance with the requirements of the regulatory region.

FEATURES

The Tilt-Pro is for use with the Portrait™ Intracranial, Head and Neck Device, AccuFix™ Head and Neck Device, and CURVE™ Device and is intended to immobilize the head and neck region. It is intended to provide optimum localization for the patient in the supine or prone positions for radiation therapy. The Tilt-Pro is intended to be used with an Aquaplast RT™ or Fibreplast® Disposable U-Frame.

OPERATING INSTRUCTIONS

INSTALLATION

INSTALLATION ON PORTRAIT & CURVE DEVICES

! NOTE ! Be sure to remove rubber feet on the bottom of the Tilt-Pro before installing onto the CURVE device.

- Firmly depress black part of locking pins into the Tilt-Pro by pressing down firmly (Fig. 1).

! NOTE ! White locking pins (P/N 2001675) are included inside the packaging and are not pre-assembled.

- Align the locking pins on the Tilt-Pro with the receiving holes on to the CURVE or Portrait Device (Fig. 2).
- Secure Tilt-Pro by pressing down on the locking pins to lock them to the device (Fig. 3). Ensure Tilt-Pro is secure to the couch top prior to use.

Fig. 1
Locking Pins


Fig. 2
Receiving Holes


Fig. 3
3 Pins

INSTALLATION

INSTALLATION ON ACCUFIX CANTILEVER DEVICE

! NOTE ! Be sure that the rubber feet are installed on the bottom of the Tilt-Pro before installing onto the Cantilever device.

! NOTE ! Fig. 4 shows an AccuFix Cantilever device prior to the installation of the Tilt-Pro.


Fig. 4

1. Align 2 protruding pins on the AccuFix Cantilever device with the 2 receiving holes on the Tilt-Pro (Fig. 5).
2. Secure Tilt-Pro by attaching and screwing down with 2 locking screws, provided with the AccuFix Cantilever device (Fig. 6). Ensure Tilt-Pro is secure to the couch prior to use.

! NOTE ! Fig. 7 is a Locking Screw (4) P/N 2001655.


Fig. 5
Protruding Pins


Fig. 6
Locking Screws

Fig. 7

SET-UP

TILT-PRO ANGLE ADJUSTMENT

The Tilt-Pro allows for adjustment from 5° to 30° in 5-degree increments.

! NOTE ! Do not adjust the height if the patient is secured in the mask.

1. Loosen knob by turning counterclockwise on each side of the Tilt-Pro to change angle adjustment (Fig. 8).

! NOTE ! Tilt-Pro upgrade contains an angle adjustment ranging from 5° to 45°. Skip to Step 3 for installation instructions of 5° to 45° angle adjustment.

2. Change angle to desired setting and secure Tilt-Pro by turning the 2 knobs clockwise and then tighten to secure (Fig. 9).
3. To remove 5° to 30° angle adjustment, remove the 3 screws at the bottom of the angle adjustment using a flat head screwdriver. This includes the black screw and the 2 white screws in Fig. 9.

! NOTE ! Do the same for the other side.

4. Unscrew by turning counterclockwise and remove the knob on each side of the Tilt-Pro (Fig. 10).
5. After all the screws and knobs have been removed, replace 30° angle adjustment with 45° angle adjustment. To attach the 45° angle adjustment, tighten the 3 screws on the bottom of each side of the Tilt-Pro as well as tighten the angle adjustment knob on both sides.


Fig. 8


Fig. 9


Fig. 10

SILVERMAN AND MOLDCARE® SUPPORT SETUP FOR SUPINE POSITIONING

1. Place a Silverman head support over the corresponding indexing pins (Fig. 11).
2. Use a MOLDCARE Cushion (if desired) in conjunction with a Silverman head support to achieve the highest level of comfort and repeatability (Fig. 12). See MOLDCARE packaging for molding instructions.

Fig. 11

Fig. 12: Formed MOLDCARE on Silverman Head Support

SET UP FOR PRONE POSITIONING

! NOTE ! Prone positioning is compatible with U-Frame masks only.

There are 2 options for prone positioning:

1. RT-4521CSA, Prone Cranial Spinal Axis Treatment Kit (Fig. 13).


Fig. 13: RT-4521CSA

! NOTE ! Refer to the Product Guide and User Manual for the Prone Cranial Spinal Axis Treatment Kit (P/N 2004216) for operating instructions, product features, specifications, warnings, cautions, and other general precautions related to the use of the Prone Cranial Spinal Axis Treatment Kit.

2. RT-4552P, Prone Cushion: place the adapter onto the two white pins.

THERMOPLASTIC SETUP

The Tilt-Pro can be used with Traditional, Slimline Disposable U-Frame Masks, or Head only Portrait S-Frame masks.

! NOTE ! Refer to the Product Guide and User Manual for the Aquaplast RT & Fibreplast Thermoplastics & IntegraBite™ (P/N 2002890) for operating instructions, product features, specifications, warnings, cautions, and other general precautions related to the use of thermoplastic masks.

U-FRAME SETUP

! NOTE ! Quick release swivel locks allow emergency patient release.

1. Place the mask over the patient, lining up the 2 pins on the Tilt-Pro with the receiving holes on the mask (Fig. 14, 15).


Fig. 14


Fig. 15

Receiving Holes

! NOTE ! Use of a new mask or previously molded mask is acceptable.

2. Secure mask using the swivel locks on the Tilt-Pro (Fig. 16).
3. Loosen swivel locks to remove the mask.



Fig. 16

Swivel Locks

S-FRAME SETUP

1. Remove swivel locks on the Tilt-Pro (Fig. 17).
2. Remove the locating pins on the Tilt-Pro using a flat head screw driver (Fig. 18).
3. Place mask over patient, lining up Portrait S-Frame pins with attachment holes on the Tilt-Pro (Fig. 19).
4. Insert pins into the Portrait S-Frame holes and press until each pin locks in place (Fig. 20).

5. To remove, pull up on each pin and lift the mask off of the patient (Fig 21).

Fig. 17


Fig. 18

Fig. 19


Fig. 20


Fig. 21

MAINTENANCE

CLEANING THE SYSTEM

The device can be cleaned with a mild, non-abrasive cleaning or disinfecting solution. To clean, apply solution to clean cloth and wipe the surface. Visually inspect the device, if it is not clean, repeat the previous cleaning steps until visually clean. Use a clean cloth moistened with water to wipe the device to remove any cleaning agent residue. To dry, wipe the device with a clean, dry cloth. The following cleaning material has been tested and found to be appropriate for cleaning the device.

- Water
- Soap and Water

DISINFECTING THE SYSTEM

To disinfect the device surface, refer to specific instructions from the cleaning agent manufacturer.

- 10% Clorox® Bleach Solution
- Isopropyl Alcohol
- Cidex® 2.4% Activated Dialdehyde Solution

Periodically check all fasteners for tightness.

Fibreplast® a registered trademark of Qfix.

AccuFix, Aquaplast RT, CURVE, IntegraBite, KVUE, Portrait, and Tilt-Pro are trademarks of Qfix.

MOLDCARE® is a registered trademark of ALCARE CO, LTD.

Cidex® is a registered trademark of Johnson & Johnson.

Clorox® is a registered trademark of The Clorox Company.

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