

# Alpine™ Stirrups

BRACHYTHERAPY WORKFLOWS WITH  
LITHOTOMY AIRSHUTTLE™ JUST GOT SMARTER

---

Alpine Stirrups are meticulously designed and feature an integrated rail clamp for easy attachment and enhanced anatomical access.

- Provide optimal pelvic access while accommodating entry into a 70cm MR bore
- Integrated rail attachment clamp ensures easy setup and prevents the loss of separate components common with standard stirrups
- Intuitive handle facilitates adjustment with minimal effort
- 20% lighter than other standard stirrups on the market\*
- Smart and comfortable boot design for easy cleaning



 **CQ MEDICAL™**  
Care Smarter

# Alpine Stirrups

## SPECIFICATIONS

- Patient Weight Limit Capacity: 500 lbs
- Stirrup Weight: 10.05 lbs
- MR Conditional up to 3T, CT compatible



RT-5100-54 Alpine Stirrups (pair)



**CONTACT US TO LEARN MORE!**



[info@CQmedical.com](mailto:info@CQmedical.com) | [CQmedical.com](https://www.CQmedical.com)

Copyright ©2024 CQ Medical. All rights reserved. Alpine and CQ Medical are trademarks of Medtec LLC.  
AirShuttle is a trademark of Qfix Systems LLC. All products may not be licensed in accordance with Canadian law.

**CQ2024P193 Rev. 1**

**\*As of 07/2024**

