

Nurofit™ Head Cushion

REDEFINING IMAGE GUIDED PRECISION



FAST, SECURE POSITIONING

The Nurofit Head Cushion ensures a secure, quick setup, saving precious time during image-guided procedures.

EXCEPTIONAL IMAGE QUALITY


The inflatable design quickly forms around the patient's head, minimizing movement and discomfort, enabling you to visualize critical details during procedures.

ENHANCED WORKFLOW

The user-friendly design, featuring an integrated inflation hand pump, means the Nurofit Head Cushion is ready to use right out of the packaging, streamlining your procedural workflow.

SEAMLESS INTEGRATION

Designed to work seamlessly with the Nurofit Premiere and most head supports, the Nurofit Head Cushion effortlessly integrates into your established workflows.



XT-4000-01C Nurofit Head Cushion



SEE THE HEAD CUSHION
IN ACTION!

Nurofit Premiere™ Head Support




SECURE AND COMFORTABLE PATIENT POSITIONING

Nurofit Premiere Head Support provides rigid, secure positioning while maintaining patient comfort with a cleanable pad. It was designed to comfortably position the patient's head during interventional neuroradiology procedures. The easily adjustable support minimizes the patient's head motion and accommodates multiple head positions.


- Compatible with a wide range of table thicknesses
- Easily attach to table
- Quick adjustment of head holder
- Adjustable median & coronal head positions
- Seamless integration of optional Nurofit Head Cushion



Compact design for improved imaging access



Adjustable for patient comfort



Compatible with table sizes
12.7 mm–44.45 mm (0.50 in–1.75 in)

XT-4000-01 Nurofit Premiere Head Support

XT-4000-01C Nurofit Head Cushion



**SEE THE HEAD SUPPORT
IN ACTION!**