

Bolus

TECH DATA SHEET



PRODUCT NAME	Bolus with Skin, Skinless Bolus
PRODUCT NUMBER	Skinless: 32052, 32102, 32152, 32202 With Skin: 31051, 31101, 31151, 31152, 31102, 31052, 31202, 31053050, 31103050, 31054040, 31104040, 31033030
GENERAL DESCRIPTION	Bolus (Bolx® & Bolx II®) is made from a solid, tissue-equivalent gel. It will not sag, deform or run. Bolus is available “skinless”, or encased in a tough layer of thin plastic “skin”.
INDICATIONS FOR USE	Designed as bolus material for use in radiation therapy, diagnostic x-ray, and ultrasound imaging procedures.
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

PRIMARY MATERIALS	Solid, homogenous, tissue-equivalent gel. Skin: Polyurethane/Synthetic Rubber		
PART NUMBER	TYPE	WEIGHT (GM)	DIMENSIONS
31033030	With Skin	278.1	30 x 30 x 0.3 cm
31053050	With Skin	772.5	30 x 50 x 0.5 cm
31103050	With Skin	1545	30 x 50 x 1.0 cm
31054040	With Skin	824	40 x 40 x 0.5 cm
31104040	With Skin	1648	40 x 40 x 1.0 cm
32052	Skinless	748.4	30 x 30 x 0.5 cm
31052	With Skin	463.5	30 x 30 x 0.5 cm
32102	Skinless	918	30 x 30 x 1.0 cm
31102	With Skin	927	30 x 30 x 1.0 cm
32152	Skinless	2177.2	30 x 30 x 1.5 cm
31152	With Skin	1390.5	30 x 30 x 1.5 cm
32202	Skinless	2812.3	30 x 30 x 2.0 cm
31202	With Skin	1854	30 x 30 x 2.0 cm
31051	With Skin	115.9	15 x 15 x 0.5 cm
31101	With Skin	231.8	15 x 15 x 1.0 cm
31151	With Skin	353	15 x 15 x 1.5 cm

PACKAGING

UNIT PACK	1 piece
UNITS / BOX	1 piece
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	May 2006