

# Choosing an All In One Solution



	Alta™ Multipurpose Device	TotalRT™ Platform
<b>DEVICE DIMENSIONS</b>	199.0 (L) x 61.0 (W) x 2.3 (H)	201.3 (L) x 61.9 (W) x 3.5 (H)
<b>DEVICE WEIGHT</b>	5.4 kg (12 lbs)	7.3 kg (16 lbs)
<b>SAFE WORKING LOAD</b>	226 kg (500 lbs)	220 kg (485 lbs)
<b>INDEXING TYPE</b>	Varian Exact   7 cm	PRODIGY 2   7 cm
<b>INFINITY EDGE™</b>	Yes	No
<b>MR SAFE</b>	Yes	Yes
<b>ELEKTA HEXAPOD COMPATIBLE</b>	Yes <sup>1</sup>	No
<b>HALCYON COMPATIBLE</b>	Yes	Yes
<b>SYMPHONY COMPATIBLE</b>	Yes (Alta AirShuttle)	No
<b>MASK ATTACHMENT STYLE</b>	Type-S	Type-S
<b>COMPATIBLE HEADRESTS</b>	Silverman <sup>2</sup> and Posifix <sup>2</sup>	Silverman and Posifix <sup>3</sup>
<b>SRS DEVICE COMPATIBLE</b>	<ul style="list-style-type: none"> <li>• SRS Adapter</li> <li>• Solstice SRS System</li> </ul>	<ul style="list-style-type: none"> <li>• SRS Adapter</li> <li>• Solstice SRS System</li> </ul>
<b>SHOULDER RETRACTOR</b>	Yes (InfinityLock Patient Handles)	Yes (UCT Patient Handles)
<b>WING BOARD COMPATIBLE</b>	<ul style="list-style-type: none"> <li>• Alta ArmShuttle Elite (Direct docking)</li> <li>• Monarch Overhead Arm Positioner <sup>4</sup></li> <li>• Pro Wing Board<sup>4</sup></li> </ul>	<ul style="list-style-type: none"> <li>• Pro Wing Board <sup>5</sup></li> <li>• Monarch Overhead Arm Positioner <sup>5</sup></li> <li>• Extended Wing Board <sup>5</sup></li> <li>• Standard Wing Board <sup>5</sup></li> </ul>
<b>VACUUM CUSHION COMPATIBLE</b>	• Vac-Lok • VacQfix	• Vac-Lok • VacQfix
<b>SBRT DEVICE COMPATIBLE</b>	<ul style="list-style-type: none"> <li>• ONEBridge</li> <li>• ONE Respiratory Belt</li> <li>• ZiFix Traverse or ZiFix</li> </ul>	<ul style="list-style-type: none"> <li>• ONEBridge</li> <li>• ONE Respiratory Belt</li> <li>• ZiFix Traverse or ZiFix</li> </ul>
<b>BREASTBOARD COMPATIBLE</b>	<ul style="list-style-type: none"> <li>• Evolve</li> <li>• Access Supine</li> </ul>	<ul style="list-style-type: none"> <li>• Evolve</li> <li>• C-Qual • C-Qual with Monarch</li> </ul>
<b>BODY THERM</b>	Yes (Direct docking)	Yes <sup>6</sup>
<b>HIP &amp; PELVIC DEVICE COMPATIBLE</b>	<ul style="list-style-type: none"> <li>• CQ Medical Bellyboards</li> <li>• FlexLock Knee &amp; Feet Support</li> <li>• KneeFix &amp; FeetFix cushions</li> </ul>	<ul style="list-style-type: none"> <li>• CQ Medical Bellyboards</li> <li>• FlexLock Knee &amp; Feet Support</li> <li>• KneeFix &amp; FeetFix cushions</li> </ul>

<sup>1</sup> HexaPOD Compatible Version RT-4800-05, <sup>2</sup> Using the Alta Silverman Adapter Plate, <sup>3</sup> Using the Headrest adapter 301035 or 301029,

<sup>4</sup> Using the Alta Locating Bar MTWBL007, <sup>5</sup> Using Type-S Lok-Bar MTWBL001, <sup>6</sup> Using the ThermFix THRMFX100