

Extended Wing Board™

TECHNICAL DATA SHEET




PRODUCT NAME	Extended Wing Board
PRODUCT NUMBERS	MTWB08T Extended Wing Board with T-Grip MTWB08U Extended Wing Board with U-Grip MTWB08 Extended Wing Board (no hand grip)
GENERAL DESCRIPTION	The Extended Wing Board is constructed of durable, lightweight ABS and provide arms-up patient comfort and support. A T-Grip or U-Grip is available for positioning preference.
INDICATIONS FOR USE	The device is intended to support and aid in positioning a patient during radiologic and other medical procedures. Devices are used to aid in the support and positioning of patients during an MRI.
MR SAFE	Yes
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS


COMPONENT	MTWB08 Extended Wing Board (no hand grip)
PRIMARY MATERIALS	Black ABS
WEIGHT (APPROX)	2.3 kg / 5.1 lbs
DIMENSIONS (CM)	61 L x 57 W x 15 H

TECHNICAL DRAWING



TECHNICAL DETAILS


COMPONENT	MTGBSUB1 T-Grip, included in MTWB08T
PRIMARY MATERIALS	Black Vinyl, White DELRIN®, White Acetal
WEIGHT (APPROX)	0.66 kg / 1.45 lbs
DIMENSIONS (CM)	34 L x 9 W x 8.5 H



TECHNICAL DETAILS

COMPONENT	MTGBSUB3 U-Grip, included in MTWB08U
PRIMARY MATERIALS	White DELRIN®, Acetal
WEIGHT (APPROX)	0.58 kg / 1.27 lbs
DIMENSIONS (CM)	14 L x 8 W x 23 H

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	Each device may have several component parts in box
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	April 2012

ACCESSORIES

SYSTEM CUSTOMIZATION	Compatible with: Silverman Headrests (MTSILVER2) Timo Headrests with Lexan Plate, Set of four (MTTIMO2L)
-------------------------	--