

New Horizon™ Prone Breastboard

TECH DATA SHEET




PRODUCT NAME	New Horizon Prone Breastboard
PRODUCT NUMBER	MTHPBB01
GENERAL DESCRIPTION	The New Horizon Prone Breastboard assists in providing comfortable, reproducible patient positioning for prone breast treatments. The modular design and variety of accessories for the New Horizon allow you to customize your patient setups, including overall height adjustments.
INDICATIONS FOR USE	The device is intended to aid in patient positioning.
MR STATUS	Not intended for use in MRI (magnetic resonance imaging).
INTEGRATION	2-Pin Indexable to couhtop
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA


TECHNICAL DETAILS

PRIMARY MATERIALS	Polyurethane, Polyethylene, Black Acetal
WEIGHT (APPROX)	27.67 kg / 61 lbs
DIMENSIONS (MM)	1483 L x 558 W x 321 H (min.)

TECHNICAL DRAWING



TECHNICAL DRAWING (BREAST CUTOUT AREA)



PACKAGING

UNIT PACK	Each
UNITS / BOX	One device - may have several component parts in box
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	June 2011

ACCESSORIES

SYSTEM CUSTOMIZATION	Patient Handles: Included with MTHPBB01 Contralateral Wedge with Movable Tray : Included with MTHPBB01 Head Cushion Covers: MTHPBB103
----------------------	---