


kVue™ Couch Top with Moving Rails Technology




TECHNICAL DATA SHEET

PRODUCT NAME	kVue Couch Top with Moving Rails Technology
PRODUCT NUMBER	RT-4551KV
GENERAL DESCRIPTION	The kVue Couchtop is designed to easily adapt to patient and treatment-specific needs through the use of interchangeable kVue Inserts.
INDICATIONS FOR USE	This device is intended to immobilize, position and reposition patients undergoing radiation therapy
MR STATUS	Not intended for use in an MR environment.
INTEGRATION	With kVue Inserts via OneTouch™ Latch
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available Upon Request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	RT-4546 kVue Couch Top with Moving Rails Technology
PRIMARY MATERIALS	Carbon Fiber, Aluminum, Steel
WEIGHT (APPROX)	56 kg / 125 lbs
MAX PATIENT WEIGHT	50 kg / 550 lbs CQ Medical Approved 200 kg / 440lbs Varian Exact™ or TruBeam™ Couch Base 155 kg / 341 lbs Varian® PRS Accessory Rail Load Rating = 45 kg / 100 lbs at 10 cm / 4 in
DIMENSIONS (CM)	See below
DOSIMETRIC PROPERTIES	Support beams are not designed to be treated through. Standard Insert = 6 mm at 6 MV Water Equivalence

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	One Couchtop
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical Website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	2010

ACCESSORIES

SYSTEM CUSTOMIZATION	<p>COMPATIBLE WITH:</p> <p>kVue Indexing Back Panel with Varian Exact compatible indexing (RT-4551KVO) Interface Plate for Elekta, required for Elekta Couch installations (RT-4551KVA-ELK) Hardware Kit for Varian TrueBeam Pedestal (RT-4551KVTBHRD) Hardware Kit for Varian Exact Couch, required for Varian Exact installations (RT-4551KVVHRD)</p> <p>KVUE INSERTS:</p> <p>kVue™ Low-e™ Insert (RT-4551KV19) kVue™ DoseMax™ Insert (RT-4551KV18) kVue™ Short Insert (RT-4551KV2) kVue™ Standard Insert (RT-4551KV13) kVue™ Universal Tip Insert (RT-4551KV8) kVue™ Encompass™ SRS Insert (RT-4600KV-01) kVue™ Portrait™ Intracranial, Head & Neck Insert (RT-4552KV) kVue™ Bos® Inserts - Standard Length (1320 mm) (RT-4535KV) kVue™ Bos® Inserts - Shorter Length (1210 mm) (RT-4535KV01) kVue™ Bos® Inserts - Longer Length (1580 mm) (RT-4535KV03) kVue™ Access 360™ G2 Package (RT-4544KV-10)</p>
----------------------	---