

kVue™ Universal Tip Insert

TECHNICAL DATA SHEET




PRODUCT NAME	kVue Universal Tip Insert
PRODUCT NUMBER	RT-4551KV8
GENERAL DESCRIPTION	kVue Inserts are customizable, treatment-specific devices designed for easy integration with the kVue Couchtop.
INDICATIONS FOR USE	This device is intended to immobilize, position and reposition patients undergoing radiation therapy.
MR STATUS	Not intended for use in an MR environment.
INTEGRATION	Insert mounts directly to the kVue couchtop through the OneTouch Latch.
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	RT-4551KV8 kVue Universal Tip Insert	
PRIMARY MATERIALS	Foam Cored Composite	
WEIGHT (APPROX)	5.9 kg / 13 lbs	
MAX PATIENT WEIGHT	249 kg / 550 lbs 38.5 kg /85 lbs @ 25 cm extended from front edge	
DIMENSIONS (MM)	1004 L x 530 W x 111 H	
DOSIMETRIC PROPERTIES	WATER EQUIVALENCE	ALUMINUM EQUIVALENCE
	6 mm @ 6 MV	~0.5 mm @ 100 kVp

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	One Device
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical Website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	2009

ACCESSORIES

SYSTEM CUSTOMIZATION	Lok-Bar (RT-4551BAR)
----------------------	----------------------