


MR Spacers & Wedges

TECHNICAL DATA SHEET




PRODUCT NAME	Spacers & Wedges
PRODUCT NUMBERS	107101MR Wedge, 2.5° 107102MR Wedge, 5° 107103MR Spacer, 20 mm 107104MR Spacer, 30 mm 107106MR Set of four includes: Wedge 2.5° and 5°, Spacer 20 and 30 mm
GENERAL DESCRIPTION	Spacers & Wedges are available in a variety of heights and angles to provide a positioning offset of the patient for unique treatment needs.
INDICATIONS FOR USE	The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.
MR SAFE	Yes
INTEGRATION & COMPATIBILITY	Attach to Posifix® baseplates, extensions or overlays. Compatible with Posifix headrests.
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

PRIMARY MATERIALS	Fiber Glass, ABS		
WEIGHT (APPROX)	0.9 kg / 2 lbs (set of 4)		
DIMENSIONS (CM)	107101MR 25 L x 12 W x 4 H 107102MR 25 L x 12 W x 6 H 107103MR 25 L x 12 W x 3 H 107104MR 25 L x 12 W x 4 H		
DOSIMETRIC PROPERTIES	BEAM DIRECTION	ATTENUATION 6 MV PHOTONS 10 X 10 FIELD SIZE (CM ²)	
	Orthogonally through the permissible radiation area of the device.	Wedge, 2.5°	2.55%
		Wedge, 2.5°	2.55%
		Spacer, 20 mm	2.64%
		Spacer, 30 mm	2.64%
	Performance Specifications are for reference only, values may vary based on setup and equipment.		

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	Each device may have several component parts in box
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	October 2024