

AccuForm™ Shim

TECHNICAL DATA SHEET



PRODUCT NAME	AccuForm Shim
PRODUCT NUMBER	MTACLSHIM
GENERAL DESCRIPTION	The Shim is an accessory for use after forming the patient's AccuForm cushion. Using this shim can counteract patient swelling and/or suspected mask shrinkage, increasing patient comfort while maintaining ideal immobilization.
INDICATIONS FOR USE	The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.
CONTRAINDICATIONS	None
MR SAFE	Yes
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	MTACLSHIM AccuForm Shim
PRIMARY MATERIALS	Cross-linked Polyolefin Foam
WEIGHT (APPROX)	0.009 kg / .02 lbs
DIMENSIONS (MM)	229 L x 93 W x 1.57 H (min.)

TECHNICAL DRAWING



PACKAGING

UNIT PACK	10
UNITS / BOX	Qty 1
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	March 2021