

# Customizable Head Support

## TECHNICAL DATA SHEET




<b>PRODUCT NAME &amp; NUMBERS</b>	<b>CHS01</b> Customizable Head Support, Posifix® <b>CHS02</b> Customizable Head Support, Type-S™
<b>GENERAL DESCRIPTION</b>	The Customizable Head Support is a shell and cushion system. The shell may be used with any customizable cushion. The unique shape of the shell was designed to encourage the cushion to return the patient to a repeatable position from fraction to fraction.
<b>INDICATIONS FOR USE</b>	The intended use of the device is to provide an additional aid to the repeat positioning of the patient for radiation or other treatment. The device forms to the back of the patient's head and eliminates the usual range of movement of non-customizable head supports.  Devices are used to aid in the support and positioning of patients during an MRI.
<b>MR SAFE</b>	Yes
<b>INTEGRATION</b>	For use with Posifix, Type-S and wide Uni-frame® baseplates
<b>FDA CLEARANCE</b>	Yes
<b>CE MARKED (EU DOC)</b>	Yes
<b>OTHER COUNTRY REGISTRATIONS</b>	Available upon request
<b>COUNTRY OF MANUFACTURE</b>	USA

## TECHNICAL DETAILS

<b>COMPONENT</b>	CHS01 Customizable Head Support, Posifix®	
<b>PRIMARY MATERIALS</b>	Thermoformed polymer, closed cell foam	
<b>WEIGHT (APPROX)</b>	0.29 kg / 0.65 lbs	
<b>DIMENSIONS (MM)</b>	267 L x 183 W x 110 H	
<b>DOSIMETRIC PROPERTIES</b> (INCLUDES CHS01 AND ACCUFORM™ MTAFL2035HF)	<b>LOCATION</b>	<b>ATTENUATION 6 MV PHOTONS 10 x 10 FIELD SIZE (cm<sup>2</sup>)</b>
	Straight up through the base (PA - posterior - anterior)	3.15%
	Side of device, 45 degree below lateral	1.28%
	Superior end, parallel to couchtop	1.68%
	Performance Specifications are for reference only, values may vary based on setup and equipment.	


## TECHNICAL DRAWING



## TECHNICAL DETAILS

<b>COMPONENT</b>	CHS02 Customizable Head Support, Type-S™	
<b>PRIMARY MATERIALS</b>	Thermoformed polymer	
<b>WEIGHT (APPROX)</b>	0.27 kg / 0.60 lbs	
<b>DIMENSIONS (MM)</b>	247 L x 183 W x 96 H	
<b>DOSIMETRIC PROPERTIES</b> (INCLUDES CHS02 AND ACCUFORM™ MTAFL2035HF)	<b>LOCATION</b>	<b>ATTENUATION 6 MV PHOTONS 10 x 10 FIELD SIZE (cm<sup>2</sup>)</b>
	Straight up through the base (PA - posterior - anterior)	2.73%
	Side of device, 45 degree below lateral	1.26%
	Superior end, parallel to couchtop	1.60%
	Performance Specifications are for reference only, values may vary based on setup and equipment.	

## TECHNICAL DRAWING



## PACKAGING

UNIT PACK	Each
UNITS / BOX	One
SHELF LIFE	n/a
WARRANTY	<a href="#">See website for Product Warranty Statement</a>
STERILIZATION	Non-Sterile
LAUNCH DATE	September 2015

## ACCESSORIES

SYSTEM CUSTOMIZATION	Any customizable cushions, including:  AccuForm™ Cushions 20 x 35 cm Chamfered (MTACL2035) 20 x 35 cm Chamfered Half Fill (MTACL2035HF)
----------------------	---