

Combifix™ 3

TECH DATA SHEET




PRODUCT NAME	Combifix 3	
PRODUCT NUMBER	151004MR	Combifix 3 for lower extremity positioning. Includes Combifix 3 Baseplate, Kneefix™ 3 Cushion, Elevation Block 3, Feetfix™ 3 Cushion and Feetfix 3 Couch Support.
	151005MR	Combifix 3 for lower extremity positioning. Includes Combifix 3 Baseplate, Low Kneefix 2 Cushion, Elevation Block 3, Feetfix 3 Cushion and Feetfix 3 Couch Support.
	302071MR	Feetfix 3 Couch Support (cushion not included)
GENERAL DESCRIPTION	The Combifix 3 is a baseplate designed to combine the Kneefix 3 and Feetfix 3 cushions, providing comfortable positioning with enhanced immobilization when treating the pelvic region and other lower extremities.	
INDICATIONS FOR USE	The device is intended for positioning the hip and lower extremities of a patient during radiotherapy and diagnostics and is used to aid in the support and positioning of patients during an MRI.	
MRI COMPATIBLE	Yes	
INTEGRATION	2 or 3-Pin indexable	
FDA CLEARANCE	Yes	
CE MARKED (EU DOC)	Yes	
OTHER COUNTRY REGISTRATIONS	Available upon request	
COUNTRY OF MANUFACTURE	USA	

TECHNICAL DETAILS

PRODUCT NAME	Combifix 3
PRODUCT NUMBER	151004MR
PRIMARY MATERIALS	PVC (Baseplate), Polyethylene (Cushions)
WEIGHT (APPROX)	20.3 lbs (9.2 kg)
DIMENSIONS (MM)	915 L x 580 W x 273 H


TECHNICAL DRAWING



TECHNICAL DETAILS

PRODUCT NAME	Combifix 3
PRODUCT NUMBER	151005MR
PRIMARY MATERIALS	PVC (Baseplate), Polyethylene (Cushions)
WEIGHT (APPROX)	20.3 lbs (9.2 kg)
DIMENSIONS (MM)	915 L x 580 W x 234 H


TECHNICAL DRAWING



TECHNICAL DETAILS

PRODUCT NAME	Feetfix Couch Support
PRODUCT NUMBER	302071MR
PRIMARY MATERIALS	PVC
WEIGHT (APPROX)	6.5 lbs (2.9 kg)
DIMENSIONS (MM)	139 L x 580 W x 139 H

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	One device, may have several component parts in box
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	Combifix 3: July 2015 Combifix 3 MR Compatible: February 2018

ACCESSORIES

SYSTEM CUSTOMIZATION	301052 Kneefix 3 301048 Low Kneefix 3 301071 Feetfix 3 301055 Elevation Blocks (Qty 1) 301056 Elevation Blocks (Qty 3) 013115 Combifix Baseplate Wall Storage 302072 Shoulder Retractors (MR Unsafe)
----------------------	---