

FlexLock™ Knee Support

TECHNICAL DATA SHEET




PRODUCT NAME	FlexLock Knee Support
PRODUCT NUMBER	RT-4800-KS
GENERAL DESCRIPTION	The FlexLock support offers customizable indexing options for positioning the lower abdomen, pelvis, and extremities.
INDICATIONS FOR USE	This device is intended to immobilize, position and reposition patients undergoing radiation therapy.
MR SAFE	Yes
INTEGRATION	Indexes to couchtop via 2-Pin Lok-Bar
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	RT-4800-KS FlexLock Knee Support
PRIMARY MATERIALS	Plastic
WEIGHT (APPROX)	3.8 kg / 8.5 lbs
MAX PATIENT WEIGHT	250 kg /550 lbs
DIMENSIONS (CM)	Small Knee Support: 45.1 L x 43.8 W x 16.1 H Large Knee Support: 45.1 L x 43.8 W x 21.1 H

TECHNICAL DRAWING



## PACKAGING

UNIT PACK	Each
UNITS / BOX	(1) Knee Support
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical Website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	2020

## ACCESSORIES

SYSTEM CUSTOMIZATION	Standard 2-Pin Lok-Bar (RT-4551BAR) Leg & Foot Strap Assembly (Set of 2) (RT-4489STRAP) Shoulder Retractor Strap (RT-4800-SR)
-------------------------	---