

MR Stretcher

MR WORKFLOWS JUST GOT SMARTER

CAPACITY & PATIENT ACCOMMODATION

- High weight capacity of 300kg / 661lbs safely accommodates a broad range of patients.
- Generous tabletop size (1930 x 595mm) and overall dimensions (up to 2010mm long, 785mm wide), providing ample space for secure patient positioning.

MRI SAFETY & COMPLIANCE


- MR Conditional up to 7 Tesla.

MOBILITY & ERGONOMIC OPERATION

- Twin-wheeled castors & adjustable push handles offer exceptional maneuverability, improving transport efficiency in tight or complex clinical environments.
- Foot-operated central and directional locking brake bar enables rapid foot-control locking/unlocking of all four wheels, including directional lock for straight-line movement control.
- Hydraulic dual-pedal height adjustment.

PRACTICAL FEATURES

- Zero Velcro on the mattress promotes improved hygiene and easier cleaning/maintenance.
- Can be paired with a Symphony® AirShuttle™ and Caddie for simplified patient positioning and transfers, even across sheet-covered surfaces.



**CONTACT US TO
LEARN MORE!**

XT-MR5502 MR Stretcher (only available in the U.S.)

info@CQmedical.com | CQmedical.com

Copyright ©2025 CQ Medical. All rights reserved. AirShuttle is a trademark of CQ Medical. CQ Medical and Symphony are registered trademarks of CQ Medical. All other trademarks are property of their respective owners. All product may not be licensed in accordance with Canadian law.
CQ2025P44 Rev. 2

 **CQ MEDICAL®**
Care Smarter