

BODY PRO-LOK ONE™ SBRT IMMOBILIZATION SYSTEM



Total Body Immobilization for Hypo-Fractionation Treatment

Body Pro-Lok ONE System provides a customizable and reproducible solution for complex SBRT patient setups.



“ We have become a referral center for SBRT patients and have employed Body Pro-Lok system for several years. It is my opinion that reproducible immobilization and abdominal compression is vital to assist in accurate targeting of lung cancers and other localized tumors. We recently upgraded our Body Pro-Lok system to the new Body Pro-Lok ONE and added a second complete system. We find it to be more flexible and efficient, especially with larger patients, without sacrificing the precision I demand.

RONALD C. MCGARRY, M.D., PH. D. UNIVERSITY OF KENTUCKY, LEXINGTON, KENTUCKY



Improved Patient Positioning and Respiratory Restriction

ONEBridge™

- Variable height adjustment with just one bridge
- Lateral offset ability allows the bridge and plate to adapt to the patient's natural posture and body type
- Tilting offset ability allows optimal patient positioning for enhanced gantry clearance
- Advanced latch engineered to release in less than three seconds for patient safety




Respiratory Belt

- Respiratory Belts provide pneumatic compression and immobilization to the abdominal region, improving accuracy and patient comfort during treatment.




Devices Easily Connect to Many Platforms


Body Pro-Lok ONE System devices are compatible with many market-leading CQ Medical platforms and couchtops.



Alta



TotalIRT



ONEPlatform

Alta™ Multipurpose Device

Alta is an all in one solution that quickly adapts positioning and immobilization solutions to your clinical setup and treatment needs, including hypofractionation.

- Features InfinityEdge™ profile and Varian Exact® Style indexing
- Combines functionality of multiple devices
- Can be extended beyond end of the couchtop for improved head and neck dosimetry
- Lightweight – 5.4 kg (12 lbs)
- CT, MR, and treatment compatible



Package Includes:		Platinum Qty	Gold Qty	Silver Qty
	Alta Device	1	1	1
	MR Locating Bars	4	4	4
	ONEBridge	2	2	1
	ONEBridge Storage Base	2	2	1
	ZiFix Traverse Abdominal System	Choose 1	1	1
	ONE Respiratory Belt	Choose 1	-	-
	Wide Respiratory Plate	1	1	1
	ComfortCare™ Cushions	3	3	3
	Hand Pump	1	1	1
	Clam-Lok Cushion	2	2	2
	Alta ArmShuttle Elite	1	1	1
	T-Shape Vac-Lok	5	-	-
	InfinityLock Hand Grip Set	1	-	-
	FlexLock Knee Support	1	1	-
	FlexLock Foot Support	1	1	-
	Kneefix™ Cushion	-	-	1
	Feetfix™ Cushion	-	-	1
	Wall Mount	1	-	-
		ALTA3000	ALTA2000	ALTA1000

TotalRT™ Platform

TotalRT is an all in one solution that quickly adapts positioning and immobilization solutions to your clinical setup and treatment needs, including hypofractionation.

- Features Universal Couchtop™ profile and Prodigy™ 2 indexing
- Large, artifact-free and homogeneous treatment area for precise dose delivery
- Combines functionality of multiple devices
- Can be extended beyond end of the couchtop for improved head and neck dosimetry
- Lightweight – 7.3 kg (16 lbs)
- CT, MR, and treatment compatible

Platinum Package (TRT3000)

FlexLock Knee & Foot Supports provide customizable indexing options

A second ONEBridge with Clam-Lok Cushion provides additional leg immobilization

ONEBridge with Wide Respiratory Plate provides improved respiratory restriction

T-Shape Vac-Lok creates a secure support around the patient

MR Lok-Bars™ allow repeatable indexing of devices

Type-S™ Lok-Bar allows for Monarch indexing

Monarch™ Overhead Arm Positioner provides overhead arm positioning and allows for Vac-Lok reproducibility

ZiFix Traverse™ Abdominal System and ONE Respiratory Belt provide pneumatic compression to the abdominal region

Patient Handles provide alternative arm positioning

Package Includes:			Platinum Qty	Gold Qty	Silver Qty
	TotalRT Platform		1	1	1
	MR Lok-Bar		4	4	4
	ONEBridge		2	2	1
	ONEBridge Storage Base		2	2	1
	ZiFix Traverse Abdominal System		Choose 1	1	1
	ONE Respiratory Belt		Choose 1	-	-
	Wide Respiratory Plate		1	1	1
	ComfortCare Cushions		3	3	3
	Hand Pump		1	1	1
	Clam-Lok Cushion		2	2	2
	Monarch Overhead Arm Positioner		1	-	-
	Type-S Lok-Bar		1	1	1
	T-Shape Vac-Lok		5	-	-
	Patient Handle Set		1	-	-
	FlexLock Knee Support		1	1	-
	FlexLock Foot Support		1	1	-
	Kneefix Cushion		-	-	1
	Feetfix Cushion		-	-	1
	Overlay Stand		1	-	-
			TRT3000	TRT2000	TRT1000

ONEPlatform™

ONEPlatform provides a simple platform solution for complex SBRT patient setups.

- Features Universal Couchtop profile and Prodigy 2 indexing
- Lightweight – 7.5 kg (16.5 lbs)
- CT, MR, and treatment compatible

Platinum Package
(ONE3000)

FlexLock Knee & Foot Supports provide customizable indexing options

A second **ONEBridge with Clam-Lok Cushion** provides additional leg immobilization

ONEBridge with Wide Respiratory Plate provides improved respiratory restriction

T-Shape Vac-Lok creates a secure support around the patient

MR Lok-Bars allow repeatable indexing of devices

Monarch Overhead Arm Positioner provides overhead arm positioning and allows for Vac-Lok reproducibility

ZiFix Traverse™ Abdominal System and ONE Respiratory Belt provide pneumatic compression to the abdominal region

Patient Handles provide alternative arm positioning

Package Includes:		Platinum Qty	Gold Qty	Silver Qty
	ONEPlatform	1	1	1
	MR Lok-Bar	4	4	4
	ONEBridge	2	2	1
	ONEBridge Storage Base	2	2	1
	ZiFix Traverse Abdominal System	Choose 1	1	1
	ONE Respiratory Belt	Choose 1	-	-
	Wide Respiratory Plate	1	1	1
	ComfortCare Cushions	3	3	3
	Hand Pump	1	1	1
	Clam-Lok Cushion	2	2	2
	Monarch Overhead Arm Positioner	1	-	-
	T-Shape Vac-Lok	5	-	-
	Patient Handle Set	1	-	-
	FlexLock Knee Support	1	1	-
	FlexLock Foot Support	1	1	-
	Kneefix Cushion	-	-	1
	Feetfix Cushion	-	-	1
	Overlay Stand	1	-	-
		ONE3000	ONE2000	ONE1000

Couchtop Compatibility

Universal Couchtop™

The Universal Couchtop is a modern and comprehensive solution for enabling advanced treatments in your clinic. Specifically designed to answer the demands of IGRT delivery, the Universal Couchtop Long Extension contains no rails or junctions in the treatment area to help facilitate better dose delivery and target coverage.


The Universal Couchtop is available for LINAC and CT, providing an environment in simulation that is reproducible and dosimetrically matched in treatment.



kVue™ One Couchtop Low-e™

The kVue One Couchtop is a one-piece modular design with a homogeneous composite base and no moving rails. It offers the ability to interchange treatment specific kVue inserts that mount directly into the top surface of the couch through the One Touch Latch. The light weight kVue One is compatible with versions of 6 Degrees of Freedom (DOF) linear accelerators and robotic bases.

kVue Couchtops are available for LINAC and CT providing an environment in simulation that is reproducible and provides for more accurate dose planning and delivery.




Elekta iBEAM® evo Couchtop

- ONEBridge and ONE Respiratory Belt are compatible via a special adapter
- Device packages are available, contact your sales rep for details



Schedule Your CE Webinar or Training Course Today!



The CQ Medical Education Center is designed to **simulate a clinical environment** and **facilitate virtual webinars**. The center allows us to provide **uninterrupted support** to our clinical partners around the world.



**VIEW OUR CE WEBINAR OPPORTUNITIES OR
SCHEDULE A TRAINING COURSE**



info@CQmedical.com | CQmedical.com



Copyright ©2025 CQ Medical. All rights reserved. Body Pro-Lok ONE, Clam-Lok, ComfortCare, CQ Medical, Feetfix, Kneefix, Lok-Bar, Monarch, ONEBridge, ONEPlatform, Prodigy, TotalRT, Universal Couchtop, Vac-Lok, and Wing Board are trademarks of Medtec LLC. Alta, ArmShuttle, DoseMax, FlexLock, InfinityEdge, InfinityLock, kVue, Low-e, and ZiFix Traverse are trademarks of Qfix Systems LLC. All other trademarks are property of their respective owners.

This brochure contains information about products which may or may not be available in particular countries. These products may not be licensed in accordance with Canadian law. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides.