

Body Pro-Lok ONE™ Rails-Only for Varian Exact® and IGRT Style Couchtops




TECHNICAL DATA SHEET

PRODUCT NAME	Body Pro-Lok ONE Varian Exact/IGRT Rails-Only
PRODUCT NUMBER	BPL020
KITS THAT INCLUDE THIS PRODUCT	BPL300
GENERAL DESCRIPTION	Allows for connection of Body Pro-Lok ONEBridge™, ONE Respiratory Belt and other Body Pro-Lok ONE devices to aid in providing comfortable, total body immobilization for hypo-fractionated treatments
INTENDED USE	The intended use of this device is to support and aid in positioning a patient during radiologic and other medical procedures
MR STATUS	Not intended for use in an MR environment
INTEGRATION	Clamps to the IGRT couchtop profile and Exact Indexing CT overlays with an edge thickness of 0.76”-0.975” (19.3 mm-24.8 mm)
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	BPL020 Body Pro-Lok ONE Varian Exact/IRGT Rails-Only
PRIMARY MATERIALS	Carbon Fiber, Black Acetal, Aluminum, Stainless Steel
WEIGHT (APPROX)	2.49 kg (5.5 lbs)
DIMENSIONS (CM)	36.7 L x 61.7 W x 8.1 H
ATTENUATION	Carbon Fiber Crossbar: 6 MV: 4.6% 18 MV: 30%

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	One device – may have several component parts in box
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical website for Product Warranty Statement
STERILIZATION	Non-Sterile
RELEASE DATE	November 2016

COMPATIBLE ACCESSORIES

SYSTEM CUSTOMIZATION	ONEBridge: BPL001 ONE Respiratory Belt: BPL026
-------------------------	---