

Respiratory Plates

TECHNICAL DATA SHEET




PRODUCT NAME	Wide Respiratory Plate	Standard Respiratory Plate
PRODUCT NUMBER	BPL008	MTSBRT003
GENERAL DESCRIPTION	When attached to the ONEBridge™, a respiratory plate places pressure at the level of the diaphragm to assist in restricting movement that influences tumor movement.	
INTENDED USE	BPL008: The device is indicated to assist in the proper positioning of patients for radiation therapy including electron, photon and proton treatments, and simulation including CT and MR image acquisition. MTSBRT003: The device is intended to aid in support and positioning of patients during MR, radiological and other procedures.	
MR SAFE	Yes	
INTEGRATION	Body Pro-Lok ONEBridge	
FDA CLEARANCE	Yes	
CE MARKED (EU DOC)	Yes	
OTHER COUNTRY REGISTRATIONS	Available upon request	
COUNTRY OF MANUFACTURE	USA	


TECHNICAL DETAILS

COMPONENT	Wide Respiratory Plate	Standard Respiratory Plate
PRIMARY MATERIALS	Black Acetal, Black G10	Black Acetal, Black G10
WEIGHT (APPROX)	0.7 lbs (0.32kg)	0.8 lbs (0.36kg)
DIMENSIONS (CM)	25.8 L x 14.4 W x 13.6 H (max)	17.1 L x 13.2 W x 13.0 H (max)

TECHNICAL DRAWING - WIDE RESPIRATORY PLATE



TECHNICAL DRAWING - STANDARD RESPIRATORY PLATE



PACKAGING

COMPONENT	Wide Respiratory Plate	Standard Respiratory Plate
UNIT PACK	1 Piece	1 Piece
UNITS / BOX	1 Piece	1 Piece
SHELF LIFE	Not Applicable	Not Applicable
WARRANTY	See CQ Medical website for Product Warranty Statement	
STERILIZATION	Non-Sterile	Non-Sterile
RELEASE DATE	December 2016	2008

ACCESSORIES

COMPONENT	Wide Respiratory Plate	Standard Respiratory Plate
SYSTEM CUSTOMIZATION	ComfortCare™ Cushion BPL007	Not Applicable