


kVue™ DoseMax™ Insert

TECHNICAL DATA SHEET




PRODUCT NAME	kVue DoseMax Insert
PRODUCT NUMBER	RT-4551KV18
GENERAL DESCRIPTION	Low electron generating insert combines the rigidity of a solid panel insert with skin-sparing properties comparable to a grid panel insert.
INDICATIONS FOR USE	This device is intended to immobilize, position and reposition patients undergoing radiation therapy
MR STATUS	Not intended for use in an MR environment.
INTEGRATION	Insert mounts directly to the kVue couchtop through the OneTouch Latch.
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	RT-455KV18 kVue DoseMax Insert
PRIMARY MATERIALS	Cored Composite
WEIGHT (APPROX)	4 kg / 9 lbs
MAX PATIENT WEIGHT	249 kg / 550 lbs
DIMENSIONS (MM)	1325 L x 530 W x 28 H
DOSIMETRIC PROPERTIES	Water Equivalence: 2.0 mm

TECHNICAL DRAWING



PACKAGING

UNIT PACK	Each
UNITS / BOX	One Device
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical Website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	2015

ACCESSORIES

SYSTEM CUSTOMIZATION	COMPATIBLE WITH: Replacement Film Skin for kVue DoseMax Insert - 3 pack (RT-4551KV5F)
-------------------------	---