


COMING SOON

MR Stretcher with AirDrive™ Technology

PATIENT HANDLING
JUST GOT SMARTER



 **CQ MEDICAL**®
Care Smarter

ADVANCING PATIENT CARE

Compassionate care requires equipment that adapts to your patients' needs. The MR Stretcher is powered by a compact AirDrive technology system to provide a seamless method for patient transfers. The system accommodates a broad range of patients with a weight capacity of 227kg (500lbs)*.

With three flexible air-source control options—including a remote hand controller and intuitive foot operation—you can focus on what matters most: your patient.



Seamless Mobility and Thoughtful Design

EFFICIENT MANEUVERABILITY

Twin-wheeled castors and adjustable push handles offer smooth, controlled steering.

SECURE BRAKING

A foot-operated central and directional locking brake bar enables rapid locking and unlocking of all four wheels, including a directional lock for straight-line movement control.

EASY ADJUSTMENTS

Hydraulic dual-pedal height adjustment helps protect clinicians from strain.

PRACTICAL ACCESSORIES

An IV Pole is included, alongside an oxygen cylinder holder** integrated into the base of the stretcher for secure storage.

INFECTION CONTROL

The mattress' zero-Velcro design allows for easier cleaning between patients

Versatile Imaging Compatibility

Designed for use in a wide array of imaging modalities, including:

CT | X-Ray | Ultrasound | Fluoroscopy | MRI (MR Conditional up to 3 Tesla)



*Weight capacity of Lithotomy AirShuttle is 204kg (450lbs).

**Warning: Operator must ensure the oxygen cylinder is safe to use in your MR environment.

DESIGNED TO POWER SYMPHONY® AIRSHUTTLES™

The MR Stetcher with AirDrive Technology is the foundation of your transfer system. Select a specialized Symphony AirShuttle patient platform that aligns with your imaging, treatment and procedural environments.



Lithotomy AirShuttle

Alta AirShuttle

Lumen AirShuttle

Iris AirShuttle



Shown with Alpine™ Stirrups (RT-5100-54).
Other compatible stirrups are available.

MR Stetcher

XT-5100-ZW MR Stetcher with AirDrive Technology

XT-MR5502 MR Stetcher (only available in the U.S.)

MR Stetcher with AirDrive Technology Packages

RT-5100-500ZW-A/B/D w/ Lithotomy AirShuttle

RT-5100-700ZW w/ Alta AirShuttle

XT-5100-800ZW w/ Lumen AirShuttle

XT-5100-600ZW w/ Iris AirShuttle



info@CQmedical.com | CQmedical.com



Copyright ©2026 CQ Medical. All rights reserved. AirShuttle and Alpine are trademarks of CQ Medical. CQ Medical and Symphony are registered trademarks of CQ Medical. All other trademarks are property of their respective owners. All product may not be licensed in accordance with Canadian law.

CQ2025P44 Rev. 3