


NEW

SGRT-Optimized, Open Neck Thermoplastic Mask Designed for Tracheostomy Patients



OPEN VIEW, OPEN NECK DESIGN

Optimized for compatibility with Surface Guided Radiation Therapy (SGRT) systems. Improves patient accessibility and comfort, particularly for patients with tracheostomies, neck swelling, or concerns related to mask confinement.

FIBREPLAST® & NON-STICK THERMOPLASTIC

Designed with greater rigidity and reduced shrinkage than standard thermoplastics to support enhanced immobilization stability and reproducible patient positioning throughout treatment.

INTEGRATED SHIM™ TECHNOLOGY

Enables patient-specific customization through precise adjustment of individual mask sections for an optimized fit and comfortable immobilization.

VARIABLE PERF MATERIAL

Strategically designed perforation pattern balances rigidity, visibility, and patient comfort.

RT-A182KYSDN Type-S Assure Open View, Open Neck, Integrated Shim, Fibreplast, Variable Perf, 3.2 mm, Qty 1*



*Must purchase in multiples of 5.

info@CQmedical.com | [CQmedical.com](https://www.CQmedical.com)

Copyright ©2026 CQ Medical. All rights reserved. Integrated Shim is a trademark of CQ Medical. CQ Medical and Fibreplast are registered trademarks of CQ Medical. All products may not be licensed in accordance with Canadian law.
PN-ASSURE-OPEN-NECK

 **CQ MEDICAL®**
Care Smarter