


Precision You Can Prove with Integrated Shim™

UPGRADE THE MASK.
TRANSFORM THE OUTCOME.



🎯 ENHANCED SETUP ACCURACY

The Integrated Shim design **reduces setup errors by an average of 0.7 mm laterally and 0.32 degrees in the X-axis, outperforms standard mask pin systems***, and supports tighter treatment margins and enhanced patient safety.

🕒 SAVE TIME, TREAT MORE PATIENTS

With an **average time savings of 10 minutes per patient**, using an Integrated Shim mask streamlines setup and enhances clinical workflow efficiency.

😊 IMPROVED PATIENT EXPERIENCE

By reducing treatment time, using an Integrated Shim mask can **alleviate patient fear**, improve compliance, and minimize the need for mid-treatment repositioning — especially beneficial for long or multi-fraction sessions.

💰 OPERATIONAL EFFICIENCY & ROI

By reducing setup errors, repeat imaging, and re-simulations, the Integrated Shim mask **helps cut operational costs, reduces machine wear, and boosts overall productivity**, maximizing your return on investment.

🔄 MULTI-DIRECTIONAL PRECISION

Integrated Shim masks demonstrate **consistently high accuracy** in the longitudinal, vertical, and X-axis directions, ensuring precise patient alignment throughout treatment.

⚙️ SEAMLESS INTEGRATION


Lower initial setup errors with an Integrated Shim mask **reduces reliance on 6DOF corrections**, allowing seamless integration with IGRT, CBCT, and adaptive therapy protocols.

*COMPARED TO THE STANDARD TYPE-S™ PIN SYSTEM

Type-S™ Integrated Shim™ Masks

TO ADJUST SHIM HEIGHT


- Pull up to unlock
- Rotate to the desired height adjustment
- Depress to lock



2 mm



0 mm




4 mm

OUR BEST THERMOPLASTIC MATERIALS

Fibreplast & Aquaplast RT with Integrated Shim


VARIABLE PERF

RT-A889KYSD	Fibreplast, Head Only, 3.2 mm
RT-A892KYSD	Fibreplast, Head & Shoulder, 3.2 mm
RT-A892YSD	Aquaplast RT, Head & Shoulder, 3.2 mm



STANDARD PERF

RT-A782KSD	Fibreplast, Head & Shoulder, 2.4 mm
------------	-------------------------------------



ASSURE OPEN VIEW, MICRO PERF


RT-A876KSDGLF	Fibreplast, Head Only, Split frame, 3.2 mm
---------------	--



visionrt | COMPATIBLE

ASSURE OPEN VIEW, VARIABLE PERF

RT-A892KYSDF	Fibreplast, Head & Shoulder, 3.2 mm
--------------	-------------------------------------



visionrt | COMPATIBLE

COMPATIBLE WITH THE TYPE-S™ HEAD & NECK POSITIONING SYSTEM
ALL THERMOPLASTICS ARE MR SAFE

info@CQmedical.com | CQmedical.com

Copyright ©2025 CQ Medical. All rights reserved. Integrated Shim and Type-S are trademarks of CQ Medical. CQ Medical is a registered trademark of CQ Medical. All products may not be licensed in accordance with Canadian law.
CQ2025P38 Rev. 0