


# VacQfix™ Cushions

## TECH DATA SHEET











PRODUCT NAME	VacQfix Cushions
PRODUCT NUMBER	See Technical Details on next page
GENERAL DESCRIPTION	Custom formed for patient comfort and positioning accuracy, becoming rigid when vacuum is applied. They can be cleaned and reused after each course of treatment. Our fine denier fabric and small-diameter foam bead fill combine to increase conformity and reproducibility.
INDICATIONS FOR USE	The device is indicated to aid in supporting and positioning patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.
MR STATUS	MR Conditional – Product poses no known hazard in 1.5T and 3.0T MR environments. VacQfix valve creates minimal image artifact.
INTEGRATION	VacQfix Indexer, 2 and 3-Pin Indexable to Couchtop (cushions with indexing bars)
APPROVED CLEANING SOLUTIONS	Sani-Cloth® AF3 Germicidal Disposable Wipe (Gray Top) Super Sani-Cloth Germicidal Disposable Wipe (Purple Top) Sani-Cloth Bleach Germicidal Disposable Wipe (Orange Top) Sani-HyPerCide Germicidal Disposable Wipe (Green Top) Sani-Cloth Prime Germicidal Disposable Wipe (Pink Top)
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA












## TECHNICAL DETAILS

PRIMARY MATERIALS	Denier Fabric, Foam Beads
-------------------	---------------------------

## TECHNICAL DRAWINGS

SHAPE	CUSHION	DIMENSIONS	FILL (LITERS)	PRODUCT NUMBER	WEIGHT (APPROX)	
	Thigh Compression SBRT	40 x 80 cm	15	RT-4517-4079VF15M	2.3 lbs	1.0 kg
	Extremities	50 x 70 cm	15	RT-4517-5070F15M	2.3 lbs	1.0 kg
	Upper Body	70 x 70 cm	20	RT-4517-7070F20M	2.8 lbs	1.3 kg
	Head & Shoulders	65 x 65 cm	7	RT-4517-T6565F7M	1.4 lbs	0.6 kg
	Breast	79 x 93 cm	20	RT-4517-T7993F20M	2.8 lbs	1.3 kg
	Multipurpose	70 x 100 cm	30	RT-4517-10070F30M	3.8 lbs	1.6 kg
	Multipurpose	70 x 100 cm	35	RT-4517-10070F35M	4.4 lbs	2.0 kg
	Multipurpose	70 x 100 cm	40	RT-4517-10070F40M	4.9 lbs	2.2 kg

## TECHNICAL DRAWINGS

SHAPE	CUSHION	DIMENSIONS	FILL (LITERS)	PRODUCT NUMBER	WEIGHT (APPROX)	
	Upper Body	100 x 100 cm	30	RT-4517-100100F30M	3.8 lbs	1.6 kg
	Upper Body	100 x 100 cm	40	RT-4517-100100F40M	4.9 lbs	2.2 kg
	Trapezoid	100 x 100 cm	35	RT-4517-100100TRAPF35M	4.4 lbs	2.0 kg
	Hip or Pelvis	80 x 150 cm	65	RT-4517-80150F65M	9.6 lbs	4.3 kg
	Multipurpose	75 x 120 cm	35	RT-4517-12075F35M	4.4 lbs	2.0 kg
	Prone Breast with Headrest Cutout	55 x 96 cm	12	RT-4517-PB3M	1.9 lbs	0.9 kg
	Prone Breast without Headrest Cutout	55 x 96 cm	15	RT-4517-PB4M	2.3 lbs	1.0 kg
	Full Body	100 x 220 cm	120	RT-4517-100220F120M	13.3 lbs	6.0 kg
	Whole Body SBRT	120 tapered to 80 x 230 cm	90	RT-4517-230120T80F90M	10.1 lbs	4.6 kg
	2-Chamber	160 x 120/80 cm	90	RT-4517-0002M	10.1 lbs	4.6 kg
	3-Chamber	230 x 120/80 cm	90	RT-4517-0003M	10.1 lbs	4.6 kg

**PACKAGING**

UNIT PACK	Each
UNITS / BOX	Qty 1
SHELF LIFE	Not Applicable
WARRANTY	See CQ Medical website for Product Warranty Statement
STERILIZATION	Non-Sterile
LAUNCH DATE	