



VacQfix™ Cushions

*Custom Formed for Patient
Comfort & Positioning Accuracy*



FULL BODY



positioning
patients for life.®

www.Qfix.com

VacQfix Cushions



**Access™ Prone
VacQfix Cushion**
*See page 5 for
more information.*

Custom formed for patient comfort and positioning accuracy, becoming rigid when vacuum is applied. They can be cleaned and reused after each course of treatment. Our fine-denier fabric and small-diameter foam bead fill combine to increase conformity and reproducibility.




SingleStep™ Valve

Unlike other systems, the unique SingleStep Valve does not require a secondary cap to ensure leak-free operation.

**All VacQfix Cushions are available
with the economy 2-piece valve**

SingleStep Valve



Economy 2-Piece Valve*




SingleStep Valve is standard on all VacQfix cushions, the Economy Valve is available upon request

**(Add "M" to the end of product code when ordering).*


VacQfix Cushions Accessories



RT-4517-0003
VacQfix Cushion,
230 cm x 120 cm/80 cm,
3-Chamber



RT-4517-IND01
VacQfix Indexer



RT-4544KV-IND
Prone VacQfix
Indexer



RT-4558CB02
ZiFix™ Abdominal/Thoracic
Motion Control System




Ribs for indexing
VacQfix cushions

RT-4546C
ArmShuttle Elite

ArmShuttle™ Elite has four locating ribs which index VacQfix cushions for increased repeatability during daily setups.




VacQfix Cushions




RT-4517-4080F15 VacQfix Cushion, 40 cm x 80 cm, 15-liter fill, for extremities


*RT-4517-4079VF15 - VacQfix Cushion, Thigh Compression SBRT, 40 cm x 80 cm, Nylon, 15-liter fill is also available with the same dimensions as RT-4517-4080F15




RT-4517-5070F15 VacQfix Cushion, 50 cm x 70 cm, 15-liter fill, for extremities




RT-4517-7070F20 VacQfix Cushion, 70 cm x 70 cm, 20-liter fill, for upper body



RT-4517-T6565F7 VacQfix Cushion, 65 cm x 65 cm, T-Shaped, 7-liter fill, for head & shoulders




RT-4517-T7993F20 VacQfix Cushion, Breast Board Cushion, 79 cm x 93 cm, 20-liter fill, with SingleStep Valve



RT-4517-10070F35 VacQfix Cushion, 100 cm x 70 cm, 35-liter fill


Dimensions listed are approximate measurements for cushions with no fill. Dimensions are for reference only.

VacQfix Cushions




RT-4517-100100F40

VacQfix Cushion, 100 cm x 100 cm, 40-liter fill, for upper body




RT-4517-80150F65

VacQfix Cushion, 80 cm x 150 cm, 65-liter fill, for hip or pelvis



RT-4517-12075F35

VacQfix Cushion, 120 cm x 75 cm, 35-liter fill



RT-4517-100100TRAPF35

VacQfix Cushion, 100 cm x 100 cm, Trapezoid, 35-liter fill



RT-4544KV-IND
Prone VacQfix Indexer

RT-4517-PB3

Access Prone VacQfix Cushion with headrest cutout




RT-4544KV-IND
Prone VacQfix Indexer

RT-4517-PB4

Access Prone VacQfix Cushion without headrest cutout



Dimensions listed are approximate measurements for cushions with no fill. Dimensions are for reference only.

VacQfix Cushions




RT-4517-0002

VacQfix Cushion, 160 cm x 120 cm/80 cm,
2-Chamber



RT-4517-0003

VacQfix Cushion, 230 cm x 120 cm/80 cm,
3-Chamber





RT-4517-230120T80F90


VacQfix Cushion, Whole Body SBRT,
120 cm tapered to 80 cm x 230 cm
with SingleStep

Dimensions listed are approximate measurements for cushions with no fill. Dimensions are for reference only.

VacQfix Cushions



RT-4517-70120F40 VacQfix Cushion, 120 cm x 70 cm, 40-liter fill



RT-4517-100220F120 VacQfix Cushion, 100 cm x 220 cm, 120-liter fill, for full body

All VacQfix Cushions come with SingleStep Valve Standard

To order the Economy Valve, add "M" to the end of the corresponding product code.

RT-4517-4080F15	VacQfix Cushion, 40 cm x 80 cm, 15-liter fill, for extremities
RT-4517-5070F15	VacQfix Cushion, 50 cm x 70 cm, 15-liter fill, for extremities
RT-4517-7070F20	VacQfix Cushion, 70 cm x 70 cm, 20-liter fill, for upper body
RT-4517-T6565F7	VacQfix Cushion, 65 cm x 65 cm, T-Shaped, 7-liter fill, for head & shoulders
RT-4517-T7993F20	VacQfix Breast Board Cushion, 79 cm x 93 cm, 20-liter fill, with SingleStep Valve
RT-4517-10070F30	VacQfix Cushion, 100 cm x 70 cm, 30-liter fill
RT-4517-10070F35	VacQfix Cushion, 100 cm x 70 cm, 35-liter fill
RT-4517-10070F40	VacQfix Cushion, 100 cm x 70 cm, 40-liter fill
RT-4517-100100F30	VacQfix Cushion, 100 cm x 100 cm, 30-liter fill
RT-4517-100100F40	VacQfix Cushion, 100 cm x 100 cm, 40-liter fill, for upper body
RT-4517-80150F65	VacQfix Cushion, 80 cm x 150 cm, 65-liter fill, for hips or pelvis
RT-4517-4079VF15	VacQfix Cushion, Thigh Compression SBRT
RT-4517-70120F40	VacQfix Cushion, 120 cm x 70 cm, 40-liter fill, Gen II, SS Valve
RT-4517-12075F35	VacQfix Cushion, 120 cm x 75 cm, 35-liter fill
RT-4517-PB3	Access Prone VacQfix Cushion with headrest cutout
RT-4517-PB4	Access Prone VacQfix Cushion without headrest cutout
RT-4517-100220F120	VacQfix Cushion, 100 cm x 220 cm, 120-liter fill, for full body
RT-4517-230120T80F90	VacQfix Cushion, Whole Body SBRT, 120 cm tapered to 80cm x 230 cm with Single Step Valve
RT-4517-0002	VacQfix Cushion, 160 cm x 120 cm /80 cm, 2-Chamber
RT-4517-0003	VacQfix Cushion, 230 cm x 120 cm /80 cm, 3-Chamber

Dimensions listed are approximate measurements for cushions with no fill. Dimensions are for reference only.

Vacuum Pumps



RT-PUMP-90
VacQfix Vacuum Pump,
110-volt, Includes hose
(RT-4517-HOSE)



RT-PUMP220Q1
VacQfix Standard Dual Mode
Vacuum Pump, 220-volt
Includes hose (RT-4517-HOSE)



RT-PUMPQ-TS
Toggle Switch for
Standard Dual
Mode Pump



RT-PUMP-KIT
Hose Adapter Kit

Simple to use and allows you to mold firm, stable patient positioning cushions in minutes. The 110-Volt Dual Mode Vacuum is equipped with a foot pedal to regulate both evacuation and inflation with an audible tone, which alerts the user when the proper evacuation level is reached. The Hose Adapter Kit allows you to use your current pump with our VacQfix cushions.

RT-PUMP-90	VacQfix Vacuum Pump, 110-volt (includes hose)
RT-PUMPQ1	VacQfix Standard Dual Mode Vacuum Pump, 110-volt (includes hose)
RT-PUMP220Q1	VacQfix Standard Dual Mode Vacuum Pump, 220-volt (includes hose)
RT-4517-HOSE	VacQfix Adapter Hose (allows industry standard pumps to accept our SingleStep Valve)
RT-PUMP-KIT	Hose Adapter Kit

VacQfix Storage Rack

Lightweight, mobile and conveniently stores multiple cushion sizes. The rack is 120 cm long and includes two 30 cm extensions. The height can be adjusted between 140 cm, 152 cm, or 165 cm.

Hang VacQfix cushions on storage rack with 'S' Hooks (sold separately in packs of 10).



RT-4517-RACK	VacQfix Mobile Storage Rack - Single Bar
RT-4517-HOOKS	VacQfix 'S' Hooks to hang VacQfix Cushions on rack, pkg. of 10

2004328_EN_D

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides. © 2022 All rights reserved. Qfix is a registered trademark of Qfix. VacQfix, SingleStep, ArmShuttle, Access, and ZiFix are trademarks of Qfix.



Qfix
440 Church Road
Avondale, PA 19311
www.Qfix.com
+1 800.526.5247
sales@Qfix.com