

ADDENDUM


8002465

Loma Linda SS5 kVue™ CT Overlay

You are receiving this addendum IFU because you have ordered the Loma Linda SS5 kVue CT Overlay. Please see the notes below indicating the changes in indexing that occurs when using this device.

! NOTE ! For kVue CT Overlay with two back panels. Both back panels must be engaged with using kVue Inserts for standard indexing. Use of only the fixed back panel with kVue Inserts may result in non-standard indexing.

! NOTE ! The kVue CT Overlay weighs approximately 25 kg (55 lbs.). At least 2 people should be involved in lifting of kVue CT Overlay.



Please refer to the kVue CT Overlay IFU, P/N 2002600, for more information.



440 Church Road
Avondale, PA 19311 USA
www.Qfix.com

+1 610.268.0585 / 800.526.5247

+1 610.268.0588 / 800.831.8174

sales@Qfix.com